

Product datasheet for **MC215672**

Gm15881 (NM_001177534) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gm15881 (NM_001177534) Mouse Untagged Clone
Tag: Tag Free
Symbol: Gm15881
Synonyms: OTTMUSG00000027193
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC215672 representing NM_001177534
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGC**C

ATGGCAGACATGCCAGCAACGGTGACTCACACCACCTTCAGGAAGTCAGTATCCCCAGGCCTGGACTCTT
TGATTATGGTGGCAACCCTCGGGCACTGCTCCGGCCTCTGGCCTTCTGAGGCTACTGCAGTTGTGTT
TACCGCATGGCTTTCCCACTAGTGCTCTATATGGGTGCCTGGATGGGCCCATGGGAAGCTGGTGTGTA
TTCTCCTGGTCCCTTAGCTTCGCCATGACTGATCGTCCTGGCCGTGGAGTTGACTGGACTCCAGAGTC
GCTTGCTCCTGTGCATGGCGGGATTTCCCATCACCTGTGCCTGCTATGCTGTACTCTTGTCTGTGCATC
ATCCATCATCTATCCTATGACTTATGTGCGGTTTCATGTCTTATGATTACATCAGAAAAATATGCTGTCTGT
GCCATTATATTCTCGTGCATCTCCTGTGTGGCTTATGCGACAGAAGTGACCTGGACTCATTCAAGGCCTG
GATATACACTGGCTACATGGCTCTAAAACCTGGGATGCTCAAGGTTTTGAGTCCTTCATAGCTTGGCT
CATCTTCTCTCCTCAGCAAGCCACAACCTACGGGAAGGAACCTGCCCTGATATGGTGTCTGGTCGTG
TATATCATATGTTTTATCCTAACCATGGTGGTTATCCTGCTCTGCCTGGGAAAATGTACCAACACTTGC
CCATCCCATTACCACCTACTTGACAGGCATGACTTTCATGTGTCTCTTGTACACCACCGCCATTGT
CCTCTGGCCTCTCTACCAGTTCAATGAAGAGTTTCATGGTTTTCCCATCGAGTAATGGATTCAAGCTGC
AGTTATAAAAAGTGCCCACTCAGTGTGTCTGGGACCGCAACTGACAGTGACATTCTTGACTGGAGTCA
ACCTGATGGCATATAGCTGACTTTATGAACTCTAGCTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001177534
Insert Size: 951 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001177534.2</u> , <u>NP_001171005.1</u>
RefSeq Size:	1948 bp
RefSeq ORF:	951 bp
Locus ID:	100271928
Cytogenetics:	8 B2