

## Product datasheet for **MC215839**

### Fam71d (NM\_029069) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Fam71d (NM\_029069) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Fam71d  
**Synonyms:** 4921509E07Rik; 4930516C23Rik; GARI-L2  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC215839 representing NM\_029069  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGATGAAGAAGAACAGGAACAAGAAGGCAATGAATAAGCAAGAAGCCCTCTTTATCCCACACTATTACA  
GACTCGGGCATCTTAAGAACATATTGGATGGAGGAGAGTATGCCCTTATGTATCTCCTCTATACTGGA  
GAGCAATTTTATTCAGGTCAATAGGAGAGGTGAATCTATTTACCTGCATAACCGAGCAAAGTGGGTGACC  
GTAGGCATCTGTTCTTCCAATCCCATCTTCAAGACCCCAATATGATGCTGTTGGACATCTGACACCAG  
AAGCTCGAAAAGAACCAGAACCCCTCTTAAAAACCTTCTGAAGTCTTCCTCTTACAGGGACTCTGGTGCT  
CACCAGGTTTATCCCTCTGCAGTTTGTGACTCTCTGTGCACAGTGCTAAGAACATGCGTCTCAAGGTG  
AAGCTGATCAACGGCCGCTCCTACTACTTACAGCTCTGTGCTCCTGTGTATAAACAGGATATCATATTTT  
CTCAGTGGGTAGACCTCATCCCCCTCCTGAATCAGGAGAAAAATCAAAAATACCAAAGTCTCCGAGGTCTC  
GAGCCTTTCAGAAGTACAAAACAGCACGGACATCGCAGGCTCGATGGACATCACGGACATCACTGCTTTC  
AGAGAAGTGGCGCCAGGCACCCGAAAACACACAGATACCCTTGTAAACATCATGGAAGTGCAGACTTCT  
CAGAATATACAGATGTCAGTGTACAGATGTCAGTGTACCGGATATGCTGGCAACGGGATCAC  
AGATGTTCAAGATATAAAAATTGTACAGAAGTACAGAAGTACAGAAGTACAGAAGTACAGAAGTACAGAAGT  
TCAGAAGTCTCAGACGTCAAAATGGCCACAAACATCTCAGGGGTCAAAATAGTATTTGAAAACGATGATA  
TAATAAGACCAAGCAGGAGGAAAAGGAATATACCCTGAAACATAGGAGTCTACAAGATACAAAAGTAA  
GAGTGACTACCAAGATTCACCAATCCTGTACCATGTCAAACATAGCCCTGGCTTTGCAAGATGAAGGG  
TGTTTTCAGACTACTCAGACTCCTGTAAAAAGCAAGGAAGACATATACAAAAGAAATATGTGATGAAACTT  
CTGAAGAAAATATGATTCGTCCTTCAAACATACATTTAAAGGCCACGGAATCCAGGAGTACAAGAAGTGA  
TTCAGACATTTTCAGATGGTGAAGTGAAGCTTGAAGATATTTGCACAAGCAGTGCCTCTAGAGAACA  
TCCTTTATCCACTCTCTGGCTCACTGCCAGGGGATTTGTCTTCTTCCCTGCTGCT**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_029069
<b>Insert Size:</b>	1323 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_029069.1</a></u> , <u><a href="#">NP_083345.1</a></u>
<b>RefSeq Size:</b>	1778 bp
<b>RefSeq ORF:</b>	1323 bp
<b>Locus ID:</b>	70897
<b>Cytogenetics:</b>	12 C3