

Product datasheet for **MC216035**

Fam46a (NM_001160378) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fam46a (NM_001160378) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Fam46a
Synonyms:	BAP014; D930050G01Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >MC216035 representing NM_001160378
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCATCAGAGATACTTTTGGACTGACCAAGGCCAAGTGGCTTTCGGAGGCACTACATGGCCGAGGGCG
 AAGGGTACTTTGCCATGGCAGAGGACGAGCTGACCGCGGCCCTTACATCCCCCTGGGTGGTGACTTCGG
 CGGCGCGGCAGCAGCTTCGGCGACCGCTGCTCGGACTATTGCGAAAGCCCCACGGCGCACTGCAATGTG
 CTGAACTGGGAGCAAGTGCAGCGGTTGGACGGCATTCTGAGCGAGACCATCCCGATCCACGGGCGCGGCA
 ACTTCCCCACGCTCGAGCTGCAGCAAGCCTGATTGTGAAGGTGGTGAAGAGGCGCTGGAAGAGAAGGG
 CATAGGTGTCCGCGACGTGCGCCTCAACGGCTCAGCCGCCAGTCACGTCCTGCACCAGGACAGTGGCCTA
 GGCTACAAGGACCTGGACCTCATCTTCTGCGCTGACCTGCGTGGGAAGAGGAGTTTCAGACTGTGAAGG
 ACGTCGTTCTGGATTGCCTGTTAGACTTCTTGCCGAAGGGGTGAACAAGGAGAAGATCACACCGCTCAC
 TCTCAAGGAGGCTTATGTGCAGAAAATGGTTAAAGTGTGCAATGATTCTGACCGGTGGAGCCTTATATCC
 CTGTCAAACAACAGTGGCAAAAATGTGGAGCTGAAGTTTGTGGATTCCCTGCGGAGGCGAGTTTGAATTCA
 GTGTAGATTCTTTTCAAATCAAATTAGACTCTCTTCTCCTCTTTTACGAGTGTTACAGAGAACCCGATGAC
 TGAGACATTTACCCACCATCATCGGTGAGAGTGTCTATGGTATTTCCACGAAGCCTTTGATCACCTT
 TGCAACAAGATCATTGCCACCAGGAACCCAGAGGAAATAAGAGGGGGAGGCGCTTAAGTACTGCAATC
 TCTTAGTGCGGGGCTTCAGGCCTGCCTCCGAGGAGATCAAGACCTCCAGAGGTACATGTGTTCCCGGTT
 TTTTATTGACTTCTCAGACATCGGAGAGCAGCAGAGAAAGCTGGAGTCTACTTGAGAACCCTTCGTG
 GGCTTGAAGACCGCAAGTATGACTATCTCATGACCCTTCATGGAGTGGTAAACGAGAGTACGGTGTGCC
 TGATGGGACATGAAAGGAGGCAGACTTTAAACCTTATCACCATGCTGGCAATCCGGGTGCTAGCTGACCA
 AAATGTCATTCCCTAATGTGGCTAATGTCACTTGCTATTACCAGCCGGCCCCCTATGTTGCAGATGCTAAC
 TTTAGCAATTACTACATTGCACAAGTTCAGCCCGTGTTCACATGCCAGCAGCAGACATACTACTTGGT
 TACCCTGCAACTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja2082_f02.zip

Restriction Sites: SgfI-MluI

ACCN: NM_001160378

Insert Size: 1344 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001160378.1, NP_001153850.1</u>
RefSeq Size:	5479 bp
RefSeq ORF:	1344 bp
Locus ID:	212943
Cytogenetics:	9 E3.1