

Product datasheet for **MC216460**

Hnf4a (NM_008261) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Hnf4a (NM_008261) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Hnf4a
Synonyms:	Hnf; HNF-; HNF-4; Hnf4; Hnf4alpha; MODY; MODY1; Nr; Nr2a1; TCF-14; Tcf14
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >MC216460 representing NM_008261
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGCGACTCTCTAAAACCCTTGCCGGCATGGATATGGCCGACTACAGCGTGCCTGGACCCAGCCTACA
 CCACCCTGGAGTTTGAAAATGTGCAGGTGTTGACCATGGGCAATGACACGTCCCCATCTGAAGGTGCCAA
 CCTCAATTCATCCAACAGCCTGGCGCTAGTGCCTGTGCGCCATCTGTGGCGACCGGGCCACCGCAAA
 CACTACGGAGCCTCGAGCTGTGACGGCTGCAAGGGTTCTTCAGGAGGAGCGTGAGGAAGAACCACATGT
 ACTCCTGCAGGTTTAGCCGACAATGTGTGGTAGACAAAGATAAGAGGAACCAGTGTCTTACTGCAGGCT
 TAAGAAGTGCTCCGGGCTGGCATGAAGAAGGAAGCTGTCCAAAATGAGCGGGACCGGATCAGCACGCGG
 AGGTCAAGCTACGAGGACAGCAGCCTGCCCTCCATCAACGCGCTCCTGCAGGCAGAGTTCTGTCCCAGC
 AGATCACCTCTCCCATCTCTGGGATCAATGGCGACATTCGGGCAAAGAAGATTGCCAACATCACAGACGT
 GTGTGAGTCTATGAAGGAGCAGCTGCTGGTCTGGTTCGAGTGGGCAAGTACATCCGGCCTTCTGCGAA
 CTCTTCTGGATGACCAGGTGGCGCTGCTCAGGGCCACGCCGGTGAGCATCTGTGCTTGGAGCCACCA
 AGAGGTCCATGGTGTAAAGGACGTGCTGCTCCTAGGCAATGACTACATCGTCCCTCGGCACTGTCCAGA
 GCTAGCGGAGATGAGCCGTGTGCCATCCGCATCCTCGATGAGCTGGTCCCTGCCCTTCCAAGAGCTGCAG
 ATTGATGACAATGAATATGCCTGCCTCAAAGCCATCATCTTCTTGTATCCAGATGCCAAGGGGCTGAGTG
 ACCCGGGCAAGATCAAGCGGCTGCGGTACAGGTGCAAGTGAGCCTGGAGGATTACATCAACGACCGGCA
 GTACGACTCTCGGGCCGCTTTGGAGAGCTGCTGCTGCTGTTGCCACGCTGCAGAGCATCACCTGGCAG
 ATGATCGAACAGATCCAGTTCATCAAGCTCTTCGGCATGGCCAAGATTGACAACCTGCTGCAGGAGATGC
 TTCTCGGAGGTTCTGCCAGTGATGCACCCACACCCACCACCCCTGCACCCTCACCTGATGCAAGAACA
 CATGGGCACCAATGTATTGTTGCTAACACGATGCCCTCTCACCTCAGCAATGGACAGATGTGTGAGTGG
 CCCCAGCCAGGGGGCAGGCAGCCACTCCCGAGACTCCACAGCCATCACCACCAAGTGGCTCGGGATCTG
 AATCCTACAAGCTCCTGCCAGGAGCCATCACCACCATCGTCAAGCCTCCCTCTGCCATTCGCCAGCCAAC
 GATCACCAAGCAAGAAGCCATCTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja2614_e11.zip

Restriction Sites: SgfI-MluI

ACCN: NM_008261

Insert Size: 1425 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC039220</u> , <u>AAH39220</u>
RefSeq Size:	4371 bp
RefSeq ORF:	1425 bp
Locus ID:	15378
UniProt ID:	<u>P49698</u>
Cytogenetics:	2 84.32 cM
Gene Summary:	<p>The protein encoded by this gene is a transcription factor involved in the development of the pancreas, liver, kidney, and intestines. The encoded protein also functions to maintain glucose homeostasis. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2015]</p> <p>Transcript Variant: This variant (1) encodes the longest isoform (HNF4alpha2). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>