

Product datasheet for **MC216520**

Ttc29 (NM_183096) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	Ttc29
Synonyms:	1700031F13Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



Fully Sequenced ORF: >MC216520 representing NM_183096
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

ATGGCTACCTTTCCCCCACTTCCCATGACTCATACAAGGCTTGCGATCTTAGCCAGGCAGAAGCTACCAT
 GCTCTTCAAAAAAATCCCAAGGGCCAGTTGATCAAAGAGAAAGAGGACATAGATTATTACCTGGAACA
 GAATTTCAAAGGGCTGTCGAAAGAAGAGGTGGCTGCCACAGGAACTCATACAAGAAGCATCTGCGTG
 GACATGCTGCGAGATGGGTTCCACAAGTCTTCACCGAGCTCTTTGCACTGATGGAACAGTGGGATAAGC
 TCCGGGAAGCTGCGAAGGCACAGTCCCTCTTCTGGCAGCAGAGGCTCTGGAGGACCAGCCGACAAGCT
 AGACAATTTCTACCACTACCTGACCAGGGCCGAGGCTGCCGAGAGGAAAGGATACTATGAGGAGGTATAT
 AACAACTGTATGCTCTGGCTGTATTTCGATAATTCTGAAGACAAGTGGGTGAGGAATCACTTCTACG
 AACGGTGTTTTAATTGCTCAGTGATCAAAGCTGACGGCGGGAAGAAAGAGGCAGAGGCCGAGTCACA
 CATGGGACTCTTTTCGAGGAAGAGGTGAGCTCCTGAAGGCTGCAGAGCACTATGAAGCATTCCATGAA
 CTGACACACGACGGCTCTGGAAGGATGGTACGGGCCAGCTTCTCAACTTGGTGGCTGTGAAAGTCTCG
 TGAGGACCTACAGATTGCTCTCAGACAGGATGCTGGAAAACAAAGACTACAAACAGGCCATCAAAATCCT
 AATAAAAGCTTCTGAAATAGCCAGAGAAGGCAATGACAGAAGCATGGAAGGGGAGGCTTCTTACTACCTG
 GGCTTGGCCCACTTGGCTTCTGGAGAATATGAAACCGCACTAACAGTGCTGAATCGGTACAGTGAGATCT
 CCACAAGTCTGGATGACGATCATGGCTGGGGAGAGCCTATGAGGCTATCGCCAAGGCCCTGCAGAGCCA
 AGGGGAGACAACAGAAGCAATTAATTACCTGGAGAAGTTCGTGACAATTGCAAGAAACAACCTGCAAGC
 CTCGATATGATTCGAGCATGCACGATGCTCGGCGACATCTACAATGAAAAGGGGCAATACAGCAAGGCTT
 CCGAGTACTTCCAACAAGCTTTCAGCACCGCGATGGAGCTCATGAAACGGCGCTCATGGATGAAACCAA
 AGTCCACTATGGGATAGCCAGAGCGCACCAGATGATGCTGGCGATGAAGGGTTACATCGAATCTGCCGAC
 TCGAATGGCCTCACTGCCTGTCTGCTGGAAGGAGACTAGAACTCAAATCGAGTACGACCCGATTCTAG
 GAGAGTCTAGAAGGCCACAGAAGATAACATATATCACTCCCAGATGCTGAGGAAGAAACAAGAAGGTC
 TCCAGAAAA**CAATAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_183096

Insert Size: 1416 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_183096.3</u> , <u>NP_898919.3</u>
RefSeq Size:	2450 bp
RefSeq ORF:	1416 bp
Locus ID:	73301
Cytogenetics:	8 C1