

Product datasheet for **MC216643**

Pou2f2 (NM_011138) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pou2f2 (NM_011138) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Pou2f2
Synonyms:	Oct-2; Oct2a; Oct2b; Oct2c; Oct2d; Otf-2; Otf2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

Fully Sequenced ORF: >MC216643 representing NM_011138
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGTTCAATCCAGCATGGGGCTCCAGAATAAGAATGTCTAAGCCCCTGGAGGCCGAGAAGCAAAGTC
TGGACTCCCCGTAGAGCACACAGACACCGAAAGAAATGGACCCGACATTAACCATCAGAACCCCCAGAA
TAAAGCGTCCCCATTCTGTGTCCCAACTGGCCCCAGCACCAAGATCAAGGCTGAAGACCCCACTGAGC
GATTCAGCCCCAGCAGCACCCCCGCCCCAGCCGGCTCAGCCTCATCTGCCCAAGGCCCAACTCATGC
TGACGGGCAGCCAGCTAGCTGGGACATACAGCAACTCCTCCAGCTCCAGCAGCTGGTGCTTGCCCGG
CCACCACCTCCAGCCACCTGCTCAGTTCCTGCTGCCACAGGCACAGCAGAGTCAGCCAGGCCCTGCTACCA
ACGCCAAATCTATTCCAGCTACCTCAACAAACCCAGGGAGCTCTCCTGACCTCCAGCCCCGGGCTGGG
TTCCTACACAGGCTATGACTCGCCCCACGCTGCCCGACCCACACCTCTCACACCCGAGCCCCGAAATG
CTTGAGCGGCCCTCCACCCGAGGAGCCAGCGATCTGGAGGAGCTGGAACAGTTTGCTCGCACCTTC
AAGCAACGCCGATCAAGCTGGGCTTACACAGGGTATGTGGGCTGGCCATGGGCAAGCTCTATGGCA
ACGACTTCAGCCAAACGACCATTTCCCGCTTCGAGGCCCTCAACCTGAGCTTCAAGAACATGTGTAACCT
CAAGCCCCCTCTGGAGAAGTGGCTCAACGACGCAGAGACTATGTCTGTGGATTCAAGCCTACCCAGCCCCA
AACCAGCTGAGCAGCCCCAGCCTGGGTTTCGACGGGCTGCCGGGGCGGAGACGCAAGAAGAGGACCAGCA
TCGAGACGAATGTCCGCTTCGCCTTAGAGAAGAGTTTCTAGCGAACCAGAAGCCTACCTCAGAGGAGAT
CCTGCTGATCGCAGAGCAGCTGCACATGGAGAAGAGTATCCGCGTCTGGTTCTGCAACCGGCCCCAG
AAGGAGAAACGCATCAACCTTGCAGTGGGCCCCCATGCTGCCAGCCCCGGGAAAGCCGACCAGCTACA
GCCCTCACCTGGTCACACCCCAAGGGGGCGCAGGGACCTACCATTGTCCCAAGCTTCTAGCAGTCTGAG
CACAACAGTTACTACCTTATCCTCAGCTGTGGGACGCTCCATCCAGCCGGACAGCAGGAGGGGTGGG
GGTGGGGGCGGAGCTGCCCCCCCTCAATTCCATCCCCTCTGTCACTCCCCACCCCGGCCACCA
ACAGCACAAACCCGAGCCCTCAAGGCAGCCACTCGGCTATTGGCTTGTCGGGCTGAACCCAGCGCGGG
CCCTGGCCTCTGGTGGAACCTGCCCTTACCAGCCT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_011138

Insert Size: 1440 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_011138.2](#), [NP_035268.2](#)

RefSeq Size: 3297 bp

RefSeq ORF: 1440 bp

Locus ID: 18987

UniProt ID: [Q00196](#)

Cytogenetics: 7 13.73 cM

Gene Summary: Transcription factor that specifically binds to the octamer motif (5'-ATTTGCAT-3'). Regulates transcription in a number of tissues in addition to activating immunoglobulin gene expression. Modulates transcription transactivation by NR3C1, AR and PGR. Isoform OCT2.5 activates the U2 small nuclear RNA (snRNA) promoter. Isoforms OCT2.1, OCT2.2 and OCT2.3 activate octamer-containing promoters. Isoforms OCT2.4 and OCT2.5 repress some promoters and activate others. Isoform OCT2.7 is unable to bind to the octamer motif, but can still activate the beta-casein gene promoter at low levels.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) has an additional segment in the coding region, as compared to variant 1. The resulting isoform (2) is longer and has the same N- and C-termini, as compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.