

Product datasheet for **MC216860**

Hyal6 (NM_028920) Mouse Untagged Clone

Product data:

| | |
|--------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | Hyal6 (NM_028920) Mouse Untagged Clone |
| Tag: | Tag Free |
| Symbol: | Hyal6 |
| Synonyms: | 4932701A20Rik; Hyal-ps1; Hyalp1 |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Cell Selection: | Neomycin |



[View online »](#)

Fully Sequenced ORF: >MC216860 representing NM_028920
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTTTATCCAGTGGGTGACACAGCTTGGTGAAGTGGTGTCTTCGTCCTTCTGGTCGCTCCGGCAGCTC
 TGAACCCGCAATGCCTCCTGTGATCAAGGACCATCCTTTCAACGTTTTCTGGGCCGCCCTACTTTGTT
 TTGTAAAGATAATTTCAATGTAACATGAATCTTCAAGTATTTAACATCATTCCGAATCCTTTCGAACT
 CAGAGTGGATCCACTATTACCGTATTTTATCCAAAGGAACTCGGGTATTACCCCTATTTCTCTGAAGATG
 GAACATCCTTCTATGGTGAATACCTCAGAAAAGTGAACCTCTCTGAGCACCTCAGGAAAAGTGCCGGCGA
 CATTGCAGATGCTGTCACTTTGTGGAGATCGGAAGGACTTGCTGTTATAGACTGGGAAGGCTGGAGGCC
 CAGTGGGATAGGAATTGGGGCAGTAGAATGATATATAAAAACCACTCGCTAGCCTTCACTCGACACCACC
 ATCCTGACTGGGCAGAAACAAAAGTGAGAACAGCTGCCAGAAAAGATTTGAAAATGCTGGAAGAAGTTT
 TATGAATGTTACTCTCACGTTGGCTTTAGAAAATGAGACCAAAGCGTTTATGGGGCTTCTATCTATCCA
 GATTGCTACAATTATGATTACCGGATAAATCCAGAGTTCTACACAGGAAGTTGCCAGATGATGAAATTT
 TCCGCAATGACCAACTCATGTGGCTATGGGAAAAAAGTGACAGCACTTTATCCTTCCATATATTTGAGTAA
 AATATTAATAAATCAAACCTTAAATGCATTGAAATTTGTTTCCAGTGGAGAGAAGCCCTGAGGGTTGCC
 GAAATGGCTAGAAAGGACTATGTTTTGCCAGTTTTTCAATTTTTCCAGACCAATTTATTTGCAAAGTATTG
 AAGCTCTCTCAGAGGAGGATTTAGTGCACACAATTGGCGAAAAGTGTGCTGATTGGGAGCAGCAGGAATAAT
 ATTTGGGGAGGATATGAATATTCGGATACAAAGGAGACTTGCTTATCTGTGAGACAACTGTTTCATGGG
 CTATTGGTCCCTATGTTCTCAATGTAACATCGGCAGCAAAGCTCTGTAGTCAAAATCTATGTAACAGTC
 ATGGAAGATGTGTACGGAAGACACCTGAATCTTCTTTTATCTGCATATGCCTGAGGATGCCATAAGAA
 TTATGTGTCAAAGAAGGGTTTCAGATTTGTCATTCCCTCCCAAGTAAGCTGAAGACAATAATGAATATG
 AAGAATGGATTTGTGCCACTGTTACTATGGCTGGCATGGAGACTCCTGCAGATCTCACTCTCCAAATC
 TCCAAAAGAACAAAGCTCCTGCAAGTGGATTGAATTCAGCAGTTATTGTTGGCATGGCCTTATTTGTGAT
 TCTGATGAATTATTTTCTATCCCTACTACAATGGCAATTTTTCTTGAACCCTGAAGAGAAGAGAA
 ATTATCTTTCTT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_028920

Insert Size: 1485 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

| | |
|-------------------------------|---|
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_028920.2</u> , <u>NP_083196.2</u> |
| RefSeq Size: | 3452 bp |
| RefSeq ORF: | 1485 bp |
| Locus ID: | 74409 |
| Cytogenetics: | 6 A3.1 |