

Product datasheet for **MC217063**

Zscan22 (NM_001001447) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zscan22 (NM_001001447) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Zscan22
Synonyms:	AW227621; D530006B18Rik; Hkr2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >MC217063 representing NM_001001447
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGCCATCCCCAAGAGCCCTCTGAGCCAGTACTCTGGGAACAGGATAGCTTCCTTCAGGTGAAGGTGA
 AAGAGGAGGAAGAAGCTAGCGACTCCCCGAGCCAGGAATCCAGCCCTAGCTTTACTGCCACCCTGAGGC
 TGCACGATTGCGCTTCAGGCACTTCCGCTATGAGGAGGCATCCAGCCCTCATCAGGCGCTAGCCCACTG
 CGAGAGCTATGTTGCCAGTGGCTGAGGCCAGAGTCGTCTCCAAAGAACAGATGCTGGAGTTGCTGGTGC
 TGGAGCAGTTCTGGGTGCGCTGCCCTGAGATCCAGGCCTGGGTGGTCTCAGTGCCCAAAGAGTGG
 AAAGGAGGCTGCCGTTCTGGTGGAGGATATGACTCAACTGCTGGACAGGAGAGGATGGGAGCCAGGAGTG
 GAACGTGAAGAGGCAAGTTGCAAACAGAGCAATACAGATGAATTGGAGCCACCAAAAATGGCGACTGAAA
 CTGTCATGGGAAGCGTTTTGCCAAAATCTACCCTTGCCACACTTGTAACCTGAGAGCCACTCAGAGAG
 CCAGCCAGAATTCCTAGGAGCACTCTGGATGAAGTCTACTGCCCAAGAGATGGATTTAGAAAAGCTCTG
 GGGCCTCATATGGATGCCCCAAAGACCAGCCTGGCCATGAGTCTAACACCTCAGAAAACGGTTCCAACA
 TGTGGCCCAACTCCCTTCCCAAGACAAGGCTTCTCAGAGGAGAAATTTGGCCATTACTTGATAATGA
 GACAGTGCCTCCAGATACTTGCTCAGAGAAGAAATCCTCCAAGGATAGTGAATGTTTGA AAAACATCCAG
 AACACTTCAGCTCTGGAGGCCACCAGAAGAGCCACTCTCAGAAGACGCCCTATGCCTGTACTGAGTGTG
 GGAAGGTATTTCAGTCGGAGCACCCACCTTGTCCAACACCAGGTTGTTACACAGGAGCAAAGCCCATGC
 GTGTAAGAGTGTGGTAAAGCCTTCAGTCGCGTTGCCACCTGACACAGCACCTGAGGATTCATACTGGG
 GAGAAGCCTTATAATGTGAGGAATGTGGTAAAACCTTCAGCCGACGACCCACCTCAGCCAGCACCACGC
 GAGTACACACAGGCGAGCGGCCCTACGAGTGTGATACTTGC GGCAAAGCCTTCAGCCAGAGCACCACGC
 GACACAGCACAGCGCATCCACACCGGGGAGAAGCCCTACAGGTGTGATGTGTGTGGGAAAAGCCTTCAGT
 GACTGTTGAGCTCTGGTCCGGCACCTGAGAGTCCACTCTGGAGAGAAGCCTTATCAATGTAAAGACTGCC
 CGAAGGCCTTCGCACAAAGCTCCTCCCTCATTGAGCACCAGCGCACTCACACAGGTGAAAAGCCTTATAA
 ATGCACTGACTGTGGAAAGCTTTAGCCGAGCTCAGCTCTCATGGTGCACCTAAAGATCCATATCACC
 GTGACTCGAAGCACCCCTTAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

- Restriction Sites:** Sgfl-Mlul
- ACCN:** NM_001001447
- Insert Size:** 1491 bp
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001001447.3](#), [NP_001001447.1](#)

RefSeq Size: 3470 bp

RefSeq ORF: 1491 bp

Locus ID: 232878

Cytogenetics: 7 A1