

Product datasheet for **MC217591**

Espn (NM_207688) Mouse Untagged Clone

Product data:

| | |
|--------------------|---------------------------------------|
| Product Type: | Expression Plasmids |
| Product Name: | Espn (NM_207688) Mouse Untagged Clone |
| Tag: | Tag Free |
| Symbol: | Espn |
| Synonyms: | je |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Cell Selection: | Neomycin |



[View online »](#)

Fully Sequenced ORF: >MC217591 representing NM_207688
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGCAATAGCTTGAACACCGAGTCTTGTCCCGGGATCAATCCATGGACCTGGAGGCAAAAGCAGCTGG
 ACTCGGGTATGTCTCGCCAACACCACCATGTGGTCCAGCCAATGACCTTTGACCTGGGCTCGCCTAC
 TAGCACGTTCTCCAATATGACTCCTGCTCCTCCAGTCACTCCAGCAGCAAGGGGAGCGATCGAATCGA
 GGGATTCCAGGTGCAAGAGCTGCAGACTTACAGAGCTACATGGACATGCTGAACCCAGAGAAGAGCTTGC
 CTCGGGGCAAGCTAGGGAAGCCTTCCCCGCCACCCTCCACCACCACCACCACCAAGCTTCCCCGCCACC
 CCCACCACCACAGGCACCCAGCCGCCCCACCTCCACCAGGCTACCCAGCTCCCAATCCCCTGTGGGA
 CTGCATCTGAATAACATTTACATGCAGACCAAGAACAAGCTTCGCCATGTGGAGGTGGACTCGTCAAGG
 AGCCGAGCTCGGGCGACGGCCACTCGGGCTACACAGGCAGGACTCCGGGCTGCTCAGGCAGGATTCGA
 GCTGCTGCACAGGCAGGAGCTGCTCAGGCACAGCACCGGACTGCGCAGGCAGGACTCCGACCGCAAACAG
 CGCTCGTTCAGTAAACAGCCAGCACGGGGGACTACTACCGCCAGCTGGGCCGACGCCCGGGGAGCCGC
 TGGCCGCACGCCCGGCATGGCCACAGCGAGGAGGCGGCGCTGCTCCCCGGGAACACGTGCACAACCGG
 CTGCTCAGCGGACTCAAAGCGTCCAGGGAGCTGCCGCCGCCACCGCCGCCGCCGCTGCCCGAGGCC
 CTGAGTTCGCCGCCGCCGCCACCTCTGCCATCGAGGGCGGGGCGCAGCCTGCGGGCAGCGTCTGTT
 CCTCGTCTTACTGGCAAAGTGAGAGTCTGAGACACAGGAAGAGCACCAATCTTTCAACATGATGTC
 CCCAACGGGTGATAACTCAGAGCTTCTGGCTGAGATAAAGGCGGGCAAGAGCCTGAAGCCGACCCGCAG
 AGCAAGGGGCTGACAACCGTGTCTCAGGCAGTGGCAGCCAGCCTCCAGCCTGAGTACCAGCAGCCTC
 TGGTGTCACTGCGCCATCTCGGACTCGGAGCCCAACCCGCGCAGCCTCTGGGTCTCAGCCACTGCTCAA
 TGGCAGTGTGGTCCCGGACACCCTGCCACCCCGCACCTGGAGTCCATCTGGATGTGGAGGCCCTCATT
 CCCACTCTTGATGAGCAGGGCCGCCATCCCGGAGTGAAGCGCCAGGTGATGGTCCGCAAGCTGCAGC
 AGAAGATGCAGGAGGAAGAGGAGCAGCGGAGGAAGGAGGAAGAGGAGGAGGCCCGGCTCGCCAGCCTGCC
 TGCCTGGAGACGAGACATTTCTCGGAAGAAGCTGGAGGAGGAGAGGGAGCAGAAGCGAAAAGAGGAGGAG
 CGGCAAAGCTGGAGGAAATACAGAGGGCGAAAGAAGTGGAGAAAGCTGCGGACACTAGGCTACGACG
 AAGCCAAGCTCGCGCCCTGGCAGCGACAGGTCTCTTGAAGAAGGGGAGATCCCTAAG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul

ACCN: NM_207688

Insert Size: 1602 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_207688.2](#), [NP_997571.1](#)

RefSeq Size: 1890 bp

RefSeq ORF: 1602 bp

Locus ID: 56226

UniProt ID: [Q9ET47](#)

Cytogenetics: 4 82.9 cM

Gene Summary: Multifunctional actin-bundling protein. Plays a major role in regulating the organization, dimension, dynamics and signaling capacities of the actin filament-rich microvilli in the mechanosensory and chemosensory cells (PubMed:14657236, PubMed:15190118). Required for the assembly and stabilization of the stereociliary parallel actin bundles. Plays a crucial role in the formation and maintenance of inner ear hair cell stereocilia (PubMed:21455486). Involved in the elongation of actin in stereocilia (PubMed:19287378, PubMed:22264607). In extrastricular hair cells, required for targeting MYO3B to stereocilia tips, and for regulation of stereocilia diameter and staircase formation (PubMed:26926603).[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (2) has multiple differences in the coding region but maintains the reading frame, compared to variant 1. This variant encodes isoform 2 which is shorter than isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data because no full-length transcript from the C57BL/6J reference strain was available. The extent of this transcript is supported by its existence in a single strain, CBA/Caj.