

Product datasheet for **MC217808**

Fibp (BC031716) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fibp (BC031716) Mouse Untagged Clone
Tag: Tag Free
Symbol: Fibp
Synonyms: 2010004G08Rik; 2010005N21Rik; 3010027N18Rik
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >BC031716
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGACCAGCGAACTAGACATTTTCGTGGGAACACGACCCTTATAGATGAAGACGTGTATCGCCTCTGGC
TGGATGGTTACTCAGTGAACGATGCAGTGGCTCTGCGAGTACGCTCCGGAATCTGGAGCAGACGGGAGC
CACCACAGGAGTCTGCAGAGCGACACCATGGACCACTACCGCACCTTTCACATGCTTGAGCGTCTGCTG
CAGCGCCGCGGAAGCTGCTGCACCAGCTCATCTTCCAGATTCCTCCCTCCCGACAGACACTCCTCATCG
AGAGGTAACACCTTTGATGAGGCCTTTGTCGGGAGGTCTTGGCAAGAAGCTGTCCAAGGGTACCAA
GAAAGACCTGGATGACATCAGCACAAAACAGGAATTACTCTCAAGAGCTGCCGGAGGCAGTTTGACAAC
TTTAAACGAGTCTTCAAGGTGGTGAAGAAATGCGGGGCTCCCTGGTGGACAACATCCAGCAGCACTTCC
TCCTCTGACCGGCTAGCCAGGGATTACGCAGCCATCGTCTTTTTTCCCAACAACCGCTTTGAAACAGG
AAAGAAAAGCTGCAGTACCTGAGCTTTGGGGACTTTGCCTTCTGTGCAGAGCTTATGATCCAGAAGTGG
ACCCTTGGAGCCGTCGGTGAAGCCCACTGACCCAGACTCTCAGGTGGACGACATGGATGTGGACTTAG
ATAAGGAGTTTCTCAAGACTTGAAGGAGCTCAAGTTTCTCGTGGCTGACAAGGACCTCCTGGACCTGCA
TAAGAGCCTGGTGTGCACTGCCCTCCGGGAAAGCTGGGTGCTTCTCTGAGATGGAACCAACTTCAAG
AATCTGTCTCGGGGCTGGTGAACGTGGCTGCCAAGCTGACCCACAATAAGGATGTCAGAGACCTATTTG
TGGACCTCGTGGAGAAGTTGTGGAACCCTGCCGCTCTGACCACTGGCCACTGAGTGATGTGCGGCTCTT
CCTCAGCCAGTATTCAGCGTCAGTCCACTCCCTGGATGGCTTCCGGCACCAAGGCACTCTGGGACCGCTAC
ATGGGCACCTCCGTGGCTGCCTTCTGCGCCTCTATCATGATTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: BC031716



[View online »](#)

Insert Size:	1095 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC031716</u> , <u>AAH31716</u>
RefSeq Size:	1234 bp
RefSeq ORF:	1094 bp
Locus ID:	58249
Cytogenetics:	19 A
Gene Summary:	May be involved in mitogenic function of FGF1 (By similarity). May mediate with IER2 FGF-signaling in the establishment of laterality in the embryo (By similarity).[UniProtKB/Swiss-Prot Function]