

Product datasheet for **MC217815**

Cldn12 (NM_001193661) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cldn12 (NM_001193661) Mouse Untagged Clone
Tag: Tag Free
Symbol: Cldn12
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC217815 representing NM_001193661
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGCTGCCGAGATGTCCACGCAGCCACCGTCCTGTCTTCTGTGTGGTATTGCCTCTGTGCGCAGGCC
TCTTTGCGGGACTCTGCTTCTAACTGGAGGAAACTGCGGCTGATCACATTCAACAGAAACGAGAAGAA
CCTGACGATTTACACGGCCTGTGGTGAAGTGTGCCCGGTATGATGGAAGCAGTGACTGCCTGATGTAC
GACCGTACGTGGTACCTGTCGGTTGACCAGCTGGACCTGCGTGTCTCCAGTTTGCCTGCCTCTCAGCA
TCGTGATCGCAATGGGTGCCTTGTACTCTGCCTGATTGGAATGTGTAACACGGCCTTCAATTCTCCGT
GCCTAACATCAAAGTGGCCAAGTGTCTGGTCAATAGTGCAGGCTGCCACCTGGTGGCCGACTCCTGTTT
TTTCTGGCAGGTACCGTGAGCCTCTCTCCGTCCATCTGGGCCATCTTTATAACAGCCATCTCAACAGGA
AGTTTGAGCCGGTCTTTACCTTTGACTATGCAGTATTTGTCACTATTGCTAGCTCAGGGGTCTGTTTAT
GACTGCTCTCTGCTGTTTCGTTTGGTATTGTGCATGCAAGTCTTTGTCTCTCTTTCTGGCAACCGCTG
TACTCTACGCTCCCGGATGCACACTTACTCACAGCCCTATTCATCACGGTCCCGCCTCTCTGCCATTG
AAATCGACATTCCAGTAGTCTCACAGCACT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
ACCN: NM_001193661
Insert Size: 735 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001193661.1 , NP_001180590.1
RefSeq Size:	3639 bp
RefSeq ORF:	735 bp
Locus ID:	64945
UniProt ID:	Q9ET43
Cytogenetics:	5 A1
Gene Summary:	<p>This gene encodes a member of the claudin family. Claudins are integral membrane proteins and components of tight junction strands. Tight junction strands serve as a physical barrier to prevent solutes and water from passing freely through the paracellular space between epithelial or endothelial cell sheets, and also play critical roles in maintaining cell polarity and signal transductions. This gene, along with several other family members, is expressed in the inner ear. The protein encoded by this gene and another family member, claudin 2, are critical for vitamin D-dependent Ca²⁺ absorption between enterocytes. Multiple alternatively spliced transcript variants encoding the same protein have been found. [provided by RefSeq, Oct 2011]</p> <p>Transcript Variant: This variant (4) differs in the 5' UTR, as compared to variant 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>