

Product datasheet for MC217836

Ceacam12 (NM_001162524) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	Ceacam12
Synonyms:	1600031J20Rik; Ceacam12-C1; Ceacam12-C3; Ceacam12C1; Ceacam12C3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

Fully Sequenced ORF: >MC217836 representing NM_001162524
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGATGAACCTCTGATGCACTTTCTGTAAAGACTGTACATCCTGGCAGGAGCTCCTGCTCACAGTCTCCA
TTTAAACATGCTGGCTGCTTCCCACTACTGCCAGATAAATTGAATCAGTGCCTCCCATTTGCTGTTGA
AGGGGATAATGTTCTTCTGTTGTGCAAACTTGCCAGAGAATGTTCAAACCTTTCTGGTACAAAGGA
GGTAACTGCTCAAGATGTTTGAATTTGCAAGACACGTGATAGCTACCAATTCTAGTGTGATGGGACCTG
CACACAGTGGTAGAGAGACGATGCTCAATAATGGATCTCTGATGATCAAGAATGTCACCAGAAAAGACTC
GGGATACTACACTCTACAAATACTTGATAACAACCTCAAGACGTGAAATAATGCGTGCAGAATTCTTTGTA
CAGAGACCAATTTTAGGCTTCAGAAAGCATCTACTTCTCAACTTAAATTTAGTGTGTCACCTAGGA
TTGAAGAAAATGATGATGTTCTTCTGCTGGTTTATAATCTGCCGAGAACTTCAAGGCTTTGTCTGGCA
CAAAGGGGTATTTCCAGTTGACCATTTAAGATAGCAAGCCATTCTCCTCACCAATTCACCATGCTG
GGGCGCACATATTTAGACAGACTGACCGTATGCAGTGTGGATCCCTGCTGCTCGAAGGTCTCACAGG
AAGACACAGGACTTTACAGCCTACGAACCATACCTGTGGATTTGATGTCAGAATCTGCCATTGTGTACCT
CAAAGTAAACAGTAAGTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001162524
Insert Size:	789 bp



OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001162524.1</u> , <u>NP_001155996.1</u>
RefSeq Size:	1899 bp
RefSeq ORF:	789 bp
Locus ID:	67315
Cytogenetics:	7 A2