

Product datasheet for **MC217837**

Ccdc127 (NM_001168659) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ccdc127 (NM_001168659) Mouse Untagged Clone
Tag: Tag Free
Symbol: Ccdc127
Synonyms: 0610011N22Rik; 5430420E17Rik; A430066A18; AI256405; AW050060; C130002N06
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC217837 representing NM_001168659
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAATAACTTAAACGACCCCCAAATTGGAATATCAGGCCCAATGCCAGGGCTGATGGGGTGATGGAA
GCAAGTGAATTATGCCCTATTGGTGCCAATGCTCGGATTGCCAGCATTCGTTGGATTTGGTCTAGGGA
ATCCAAAAAGAGATAGAAAAAGCCAGAAAAGCATATCATCAGAGAACAGCTGCTTCCAACAGGATCTT
GAAGCCAAATACCATGCTGTGATCTCAGAACACCGCGCGCTGTGGCTCAGCTGTCTTTGAACTGGAAA
AAGAGCAAAATAGAACAAGTAGTTTTCCGAGAAGCCCTCATCTCGCAGGGCGTAAGTTGGCAGAAGAAAA
GAACTTCTCGAACAGGAACGGGCTCAGATAAAGCAAGAAAAAGCCGGCTCCAGCCTCTGAGAAATGTC
TATTTGAGCTGCCTGCAAGAAGAGGATGACTGGCAGAGGAGAGCCAGCATGTGCTGAAAGAGGTGGGAG
AAGCTCTCGAGGAACGCCAAAATATCTATTGTAGCCTGATCATCCCCGTAGTGCACGACTAGAGCTGGA
GAAGAGCTTACTGGTGCGCACATCTGTTGATCCAGTGGCTGTGACCTAGAAATGGCAGCTGGCTTGTCT
GACATATTTAAGCATGATAAGCATTGTGGAGATGTCTGGAACACCAACAAACGTCAAAATGGGAAGCTCA
TGTGGATGTATCTCAAATACTGGAACTACTTGTGCAACTTAAAAATTTAAGAAAGTAGAGAAGGTCAT
ATTAGAAAAATAA

ACGCGTACGCGGCCGCTCGAGCAGAAAACATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001168659
Insert Size: 783 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001168659.1</u> , <u>NP_001162130.1</u>
RefSeq Size:	3436 bp
RefSeq ORF:	783 bp
Locus ID:	67433
UniProt ID:	<u>Q3TC33</u>
Cytogenetics:	13 C1