

Product datasheet for MC217925

Fam122b (NM_001166365) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	Fam122b
Synonyms:	4632404H22Rik; AI256344; SPACIA2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC217925 representing NM_001166365 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCTCAGGAGAAAATGGATCTGGATTTTGAGGCTGACACATCTGAGGGGGCCACCCTGAGACGGTCCA
ACAGTGCTCCCTGATCCATGTGCTCAGTGACCTTTCACAGGTTTTTGAGCCATATCCACTTAGAACTGG
GAGGACCAGTACGGCAATTATGAGCCATCATAGCTTGAAGAAGGCCTGGATATGATGAACAGAGAACT
ACAAATGAAAGGGAAGCACAAGCAGGGATGCAGATAAGCCAGTCATGGGACGAAAGCCTGAGCCTGAGTG
ATAGTGATTTTGACAAGCCAGAGAAATTATTTCTCCAAAAGAATTGACTTCACTCCAGTTTCTCCAGC
TCCTTCCCAACCAGAGGATTTGGAAGCAATGCTTGTGCGCCATCCTTGCAAATGTTTGTGAGCAGCAGC
GGGATGCCACCGAGTCTGTTCTCAATCCAAGACATTTTCAAGGAGGAGCCAGAGCCCGGTCAAGTGCA
TTAGGCCAGTGTTCTTGGTCTCTTAAGAGAAAAGGGAAATGGAGATGGAAGCCAGCCTAAGAGACC
CTTCAAGGCACCACCAGCATGCTGTCTACCAATCCTGCACAATTGTCTGATTTCAAGTTTCATGTTTCAGAT
ATTTTGGATGGCAGTAGTATCAGCAGTGGCTTGTCTCAGACTCACTGGCTACAGGCAGCGTCCCGCAG
AGTCTCCAGTAGCATGCTCCAATTCATGCTCTCCGTTTCATCTTGATGGATGATCTCTACCCAAGTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001166365
Insert Size:	768 bp



OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001166365.2</u> , <u>NP_001159837.2</u>
RefSeq Size:	4331 bp
RefSeq ORF:	768 bp
Locus ID:	78755
UniProt ID:	<u>Q6NZE7</u>
Cytogenetics:	X A5