

Product datasheet for MC217964

Hexim1 (BC022111) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Hexim1 (BC022111) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Hexim1
Synonyms:	HIS1, CLP-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC217964 representing BC022111 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCAGAGTCCCTCGTGCAGCCAGGTCATGACTCGGAGGCCACCAAGCAGGAGGCTCCTGCCGCCGGAG
GCGAGGAGCCATGGGACAGCAACAGAGACAGCTGGGCAAGAAAAACATCGGAGACGCCCTCAAAGAA
GAAGCGGCATTGGAAGCCCTACTACAAGCTGACTTGGGAGGAGAAGAAAAAGTTCGACGAGAAGCAGAGC
CTGCGAGCTTCGCGGGTTCGAGCCGAGATGTTCCGAAGGGCCAGCCAGTTGCGCCCTATAACACCAGC
AGTTCCTCATGGATGACCACGATCAGGAGGAGCCTGATCTCAAAACCGGCCTTTACCCAAAGCGGGCAGC
CGCCAAATCCGACGACACCAGCGATGAGGATTTTGTGGAAGAAGCTGGTGAGGAAGCAGGAGGCAGCGAT
GGCATGGGAGGGGACGGCAGCGAGTTTCTGCAGCGGGACTTCTCGGAGACGTACGAGCGGTACCACGCCG
AGAGCCTGCAGAACATGAGCAAGCAGGAGCTCATCAAAGAGTACCTGGAGCTGGAGAAGTGCCTCTCCCG
CAAGGAGGACGAAAAAACCAGGCTGCGGCTGGAAGCAAGCGGCTGGGAGGCGTCGACGCGCGAGTGCAG
GAGCTCGAGCTAGAGCTGGACCGGCTGCGCGCTGAGAACCTCCAGCTGCTGACCGAGAACGAACTGCACC
GGCAGCAGGAGCGAGCGCCTCTTCCAAGTTCGGCGACTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	BC022111
Insert Size:	741 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC022111 , AAH22111
RefSeq Size:	1488 bp
RefSeq ORF:	740 bp
Locus ID:	192231
Cytogenetics:	11 E1
Gene Summary:	Transcriptional regulator which functions as a general RNA polymerase II transcription inhibitor. In cooperation with 7SK snRNA sequesters P-TEFb in a large inactive 7SK snRNP complex preventing RNA polymerase II phosphorylation and subsequent transcriptional elongation. May also regulate NF-kappa-B, ESR1, NR3C1 and CIITA-dependent transcriptional activity. Plays a role in the regulation of DNA virus-mediated innate immune response by assembling into the HDP-RNP complex, a complex that serves as a platform for IRF3 phosphorylation and subsequent innate immune response activation through the cGAS-STING pathway.[UniProtKB/Swiss-Prot Function]