

Product datasheet for MC218023

Fbxo27 (NM_001163702) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids

Tag: Tag Free

Symbol: Fbxo27

Synonyms: E130008B10Rik; FBG5; Gm161

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >MC218023 representing NM_001163702
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTGCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGGTGCCTGGATATCCAGGACCCGGTCCCAGCGCCGAGCCCGACCCAGGAGGTCTGGACCTGA
GCCGTCTGCCCCGGAGCTGCTCCTGCTGGTACTGAGCCACGTGCCCGCGCACGCTGCTGATGCACTG
CCGCCGAGTGTGCCGCGCTGGCGCGCCCTGGTAGATGGCCAGGCCCTATGGTTGCTGCTGCTGGCCGT
GACCACAGCGCAGCGGGCCGCGCCCTGCTGACACTAGCCCGCCGCTGTCTGCCCCCTGCCACGAAGACA
CGCCCTGCCCGCTGGGCAATTCTGCGCACTGAGACCGCTAGGACGCAACCTCATCAGCAACCCCTGTGG
CCAAGAAGCCTGCGCAAATGGATGGTGGCACGCTGGGGATGGATGGGTGGTGGAGAAGAAGAGGAAA
CCTGTACCTGGGGCCCTTCACAGACCTGCTTCGTGACTTCCTTCAGCTGGTGTGCAAGAAGCAGGTGCG
TGGACCTGGTGGAGAAGGTCTGTGGCCAGAGCTGCTGGACAGTGGTGGTGGAGATTGCTGTCTCTGA
CTGGTGGGAGCTCGACATGACAGTGGCTGTAAGTACCGTCTCTTTGTACGCTTCTTGATGCCACCAG
AACGTCATAGATAAGTTCTCCGCTGTGCCAGATCCCATTGAACAGTGAACAATGATATCTACCTGCAGG
TCACCCACGTGTTCTCCGCATCAGGAGGGGCATACGTTTTGTCTCTTTGAACACTGGGGCCAGGACAC
ACAGTTCTGGGCTGGCCACTATGGGCCAGAGTGACGAACCTCCAGTGTGATCATACGAGTCTGTCACTCC
TAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001163702



Insert Size:	843 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001163702.1, NP_001157174.1</u>
RefSeq Size:	1966 bp
RefSeq ORF:	843 bp
Locus ID:	233040
UniProt ID:	<u>Q6DIA9</u>
Cytogenetics:	7 B1
Gene Summary:	Substrate-recognition component of the SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex. Able to recognize and bind complex-type oligosaccharides.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longest transcript and encodes the longer isoform (1).