

Product datasheet for MC218058

Fhad1 (NM_028429) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fhad1 (NM_028429) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Fhad1
Synonyms:	B230311F11, 2900090M10Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC218058 representing NM_028429 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGTAGAGGACGCTCAGCAGATGAAGACAACGGAGAGTCAAAGAGCCGAGACGCTCGCTCTGAAGTTAA
AGGAAACATTAGCAGAGCTAGAAACAACCAAGACGAAAATGATCCTGACGGACGATCGGCTGAAGCTGCA
GCAGCAGTCGATGAAGGCCCTGCAGGATGAGCGGGAATCCAGAAGCACGGCTTCGAGGAGGAAATCTCA
GAGTACAAAGAGCAGATCAAACAGCACTCGCAGACCATCGTGAGCCTAGAGGAAAGGCTCTGCCAAGTGA
CTCAGTACTACCAGAAGATAGAAGGGGAGATCACAACTGTGAAAAACAACGATACAGGACCGAAGGAGGA
AGCCTCCCAAGACCTTACAGCAGGACCCCCCTGGACAGTGGTGACAAAGAAATAGCGTGCGACCATCTG
ATAGATGACTTACTAATGGCTCAGAAGGAAATCCTATCCAGCAGGAGATAATCATGAAGCTGAGGACTG
ACCTCGGGGAAGCTCACAGCCGGATGTCAGACCTGAGAGGGGAACTAAGCGAGAAGCAGAAGATGGAGCT
GGAGCGGCAGGTGGCCCTGGTGCGACAGCAGAGCGGGAGCTGAGCATGCTGAAGCGAAGGTGGCGCAG
ACGACTGGCCTGATGGAGAAGAAGGACAGGGAAGTGAAGTCTCAGGGAGGCGCTCAGGTGTGGGAAGG
GGTCCGT**GTA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_028429
Insert Size:	711 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_028429.1, NP_082705.1</u>
RefSeq Size:	882 bp
RefSeq ORF:	710 bp
Locus ID:	329977
Cytogenetics:	4 E1