

Product datasheet for MC218076

Golm1 (NM_001035122) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Golm1 (NM_001035122) Mouse Untagged Clone
Tag: Tag Free
Symbol: Golm1
Synonyms: 2310001L02Rik; AW125446; D030064E01Rik; Golph2; GP73; PSEC0257
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC218076 representing NM_001035122
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGC**C

ATGATGGGATTGGGAATGGGCGTCGCAGCATGAAGTCTCGCCCTTATCCTGGCTGCCCTGGTGGCCT
 GTGTTATTGTGCTAGGCTTCACTACTGGATTGCTAGCTCGAGAAGCGTGAACCTCAGACGCGCATTGT
 AGAGCTGGAAGGCCGGTCCGCAGGGCGGCTGCAGAGAGAGGTGCTGTAGAGCTCAAGAAGAAATGAGTTC
 CAGGGCGAGCTGCAGAAGCAGAGAGAGCAACTTGACAGGATCCAGTCTAGCCACAGCTTCCAGCTGGAGA
 ATGTCAACAAGCTACACCAGGATGAGAAGGCGGTTTTAGTGAATAACATCACCACAGGCGAGAAGCTCAT
 TCGAGACCTGCAAGACCAGTTGAAGGCACTGCAGAGGAGTTACAGCAGCCTGCAGCAGGATATCTTCCAA
 TTCCAGAAGAACCAGACCAGCCTGGAGAAGAAGTTCTCCTATGACCTGAACCAAGTGTATCAGCCAGATGA
 CCGAGGTGAAGGAACAGTGTGACGAGCGGATAGAGGAGGTCATCAGAAAGCGAAACGAAGCACCGGGTTC
 CAGAGACCTGGCTGAAACCAACAACCAACATCAGCAGGCCCTCAAGCCTCAGCCCAAGCTTCAGGAAGAA
 GTCCCATCAGAGGAACAGATGCCACAAGAGAAAGGGGATGTGCCAGAAACAAGTCCAGATTCCAGCCC
 CCAATTCAGAGAGCTTGGGTTTGAAGCCACAGGTGCAGAATGAGGAAACCAATGAGATTCCAGGCTGTCGG
 TGAGGAGCACCAACAGGCATCCATCCAAGGCGAGGCGGTGGCGGATGGCACTCGGGTAGGAGCTGAGAAA
 CTCGATCAGCACACAGCTACCAGCTGGCCTGTTGGCCAGGCCTGAGGAGGACTCTCAGTACCCTGAGC
 GGGAGCAGCTTGTATCCGGGATAGGCAGGAGCAGCAGCGTGCATCAGAGGAAGGGGAGGCCAGAAGAA
 CCCAGGAGATGAATACGACATGGATGAAAACGAAGCAGAATCTGAGAGGGAGAAGCAGGCAGCCCTTGCA
 GGGAAATGACAGAAATATAAATGTTCTTAATGCTGACGCTCAGAAGAGAGGCATAATAAATGTGCCCGTCG
 GGAGTGAAGGCAGAGCCACATTCTAAATCAAGTAGGCATTTCATATTCCTCAACAAGCC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online >](#)

ACCN:	NM_001035122
Insert Size:	1182 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001035122.2, NP_001030294.1</u>
RefSeq Size:	3557 bp
RefSeq ORF:	1182 bp
Locus ID:	105348
UniProt ID:	<u>Q91XA2</u>
Cytogenetics:	13 B2
Gene Summary:	Unknown. Cellular response protein to viral infection (By similarity).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) lacks an alternate segment in the 5' UTR, compared to variant 1. Both variants encode the same protein.