

Product datasheet for **MC218188**

Igk (BC091754) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Igk (BC091754) Mouse Untagged Clone
Tag: Tag Free
Symbol: Igk
Synonyms: 19-23; EG637227; Gm7202; IgK
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >BC091754
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGGCATCAAGATGGAGACACATTCTCAGGTCTTTGTATACATGTTGCTGTGGTTGTCTGGTATTGAAG
GAGACATTGTGATGACCCAGTCTCACAATTCATGTCCACATCAATAGGAGACAGGGTCAGCATCACCTG
CAAGGCCAGTCAGGATGTGACTACTGCTGTAGCCTGGTATCAACAGAAACCAGGACAATCTCTAAAGTA
CTGATTCAGTGGGCATCCACCCGGCACACTGGAGTCCCTGATCGCTTACAGGCAGTGGATCTGGGACAG
ATTTCACTCTCACCATTAGCAATGTGCAGTCTGAAGACTTGCCAGATTATTTCTGTGCAGCAATATAGCAG
CTATCCTCTCACGTTTCGGCTCGGGGACAAAGTTGGAAATAAACGGGCTGATGCTGCACCAACTGTATCC
ATCTTCCCACCATCCAGTGAGCAGTTAACATCTGGAGGTGCCTCAGTCGTGTGCTTCTTGAACAATTCT
ACCCAAAGACATCAATGTCAAGTGAAGATTGATGGCAGTGAACGACAAAATGGCGTCTGAACAGTTG
GACTGATCAGGACAGCAAAGACAGCACCTACAGCATGAGCAGCACCCCTCACGTTGACCAAGGACGAGTAT
GAACGACATAACAGCTATACCTGTGAGGCCACTACAAGACATCAACTTCACCCATTGTCAAGAGCTTCA
ACAGGAATGAGTGTTAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: BC091754
Insert Size: 657 bp



OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC091754
RefSeq Size:	948 bp
RefSeq ORF:	716 bp
Locus ID:	16071
Cytogenetics:	6 32.04 cM