

Product datasheet for **MC218191**

Il18r1 (NM_008365) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Il18r1 (NM_008365) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Il18r1
Synonyms:	Il1rrp; Il18ralpha
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >MC218191 representing NM_008365
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCATCATGAAGAATTAATCTTGACACTCTGCATTCTCATTGTTAAAAGTGCCTCAAAAAGTTGTATTC
 ACCGATCACAAATTCATGTGGTAGAGGGAGAACCCTTTTATCTGAAGCCATGTGGCATATCTGCACCAGT
 GCACAGGAATGAAACAGCCACCATGAGATGGTTCAAAGGCAGTGCTTCACATGAGTATAGAGAGCTGAAC
 AACAGAAGCTCGCCAGAGTCACTTTTCATGATCACACCTTGAATTCTGGCCAGTTGAGATGGAGGATG
 AGGGAACGTACATTTCTCAAGTCGGAATGATCGTCGCAATTGGACCTTAAATGTCACCAAAAAGAAACAA
 ACACAGCTGTTTCTGACAAGCTCGTGAAGCAGAGATGTTGAAGTTAACAAATCTCTGCATATCACT
 TGTAAGAATCCTAACTATGAAGAGCTGATCCAGGACACATGGCTGTATAAGAACTGAAGGAAATATCCA
 AAACCCCAAGGATCCTGAAGGATGCCGAGTTTGGAGATGAGGGCTACTACTCTCGCTGTTTTCTGTCCA
 CCATAATGGGACACGGTACAACATCACCAAGACTGTCAATATAACAGTTATTGAAGGAAGGAGTAAAGTA
 ACTCCAGCTATTTTAGGACCAAAGTGTGAGAAGGTTGGTGTAGAAGTAAAGGATGTGGAGTTGAACT
 GCAGTGCTTCATTGAATAAAGACGATCTGTTTTATTGGAGCATCAGGAAAAGGACAGCTCAGACCCATA
 TGTGCAAGAAGACAGGAAGGAGACGACAACATGGATTTCTGAAGGCAAATGCATGCTTCAAAAATACTG
 AGATTTAGAAAAATTAAGAACTATCTCAATGTTTTATATAATTGCACCGTGGCCAACGAAGAAGCCA
 TAGACACCAAGAGCTTCGTCTTGGTGAGAAAAGAAATACCTGATATCCCAGGCCATGTCTTTACAGGAGG
 AGTAACTGTGCTTGTCTCGCTCTGTGGCAGCAGTGTGTATAGTGATTTTGTGTGCATTTATAAAGTT
 GACTTGGTCTGTCTATAGGCGCATAGCGGAAAGAGACGAGACACTAACAGATGGTAAAACATATGATG
 CCTTTGTGCTTACCTGAAAGAGTGTCTCCTGAGAATAAAGAAGAGTATACTTTTGTGTGGAGACGTT
 ACCCAGGTCCTGGAGAAACAGTTTGGGTATAAGTTATGCATATTTGAAAGAGATGTGGTGCCTGGCGGA
 GCTGTTGTCGAGGAGATCCATTCATGATAGAGAAAAGCCGAGGCTAATCATCGTTCTCAGCCAGAGTT
 ACCTGACTAACGGAGCCAGGCGTGAGCTCGAGAGTGGACTCCACGAAGCACTGGTAGAGAGGAAGATTAA
 GATCATCTTAATTGAGTTTACTCCAGCCAGCAACATCACCTTTCTCCCCCGTGCCTGAAACTCCTGAAG
 TCCTACAGAGTTCTAAAATGGAGGGCTGACAGTCCCTCCATGAACTCAAGGTTCTGGAAGAATCTTGTTT
 ACCTGATGCCCGAAAAGCCGTCAAGCCATGGAGAGAGGAGTCCGAGGCGCGGTCTGTTCTCTCAGCACC
 TTA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja2147_a07.zip

Restriction Sites: SgfI-MluI

ACCN: NM_008365

Insert Size: 1614 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_008365.2](#), [NP_032391.1](#)

RefSeq Size: 4158 bp

RefSeq ORF: 1614 bp

Locus ID: 16182

UniProt ID: [Q61098](#)

Cytogenetics: 1 B

Gene Summary: Within the IL18 receptor complex, responsible for the binding of the proinflammatory cytokine IL18, but not IL1A nor IL1B. Involved in IL18-mediated IFNG synthesis from T-helper 1 (Th1) cells (By similarity). Contributes to IL18-induced cytokine production, either independently of SLC12A3, or as a complex with SLC12A3 (PubMed:26099046). [UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (a). Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments and are consistent with Northern analyses.