

Product datasheet for MC218209

Nnp1 (BC022910) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nnp1 (BC022910) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Nnp1
Synonyms:	AL033305; AL033306; AU018604; NNP-1; Nnp1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC022910 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATCCGGGAGTGGGTGGGCATTGACAGACTGCGCCTGGATAAGTTCTACATGCTCATGAGGATGGTCC
TAAGTGAGTCACTGAAGGCTGTGAAAGCCCGGGATGGGACGAGAGACAGATTGAGCAGCTGCTAGAGCT
GCTAACACCGAGATCCTGAACCCGACAGCCAGGCCCCAGTGGGTGAAGAGCCACTTCTAGAGATC
TTTCTGGAGGAGCTGGCCAAAGTGGCGCTGCAGAGCTCACCGCAGACCAGAATCTGCAGTTCATCGACC
CCTTCTGCCAGATAGCAGCCGCACTAAGGATTCCTAGGTTTGCATAAGATCATTGAGCATCTTTCA
GACGATTGTGGAGCAGGCTCCCTTGGCCATAGAAGACATCATGAATGAACTGGACACACAGAGTGGAGAG
GGTGAAGCATCAGATGGTGATGATGGTGAAGCATCAGATGGTGATGATGGAGAAGCATCAGACGATGATG
ATGGTGAAGCGTCAGATGGTGGTGATGGTGATGTGGCAGACAGTGATGATAGTGATGGGCGACAGATGA
TGATGGTGACGTGTGACACGGTGACGGTGGTGACAATGATGAAGGTGACAGTAACAAGTCTCAGAGGGT
GAGCAGGATCTGCAAGACACACCACCAAGAAGCTGCCTGCAGGCACCGCACATAGGGCTGGCCCTGAGG
CAGACAAGGAACAGGCGTGGGATGACGAGGAGAACGCTGGCCCTGTCTCCAGTTTGACTACGAGGCACT
TGCTAACAGGCTGTTCAAGCTCGCCAGTAGACAGAGCACCCCTTCTCAGAACAGGAAGCGCTCTACAAA
GTGATCCAAAAGCTGCGGAGCTGGCTGGAGGCACCTTCCAGAGGATGATGTCCCTGAAAAAGCCTACA
AGAAGATGCTGGAAGGGAGGCGGGAACGGAAGAAAGAAAAAGAGCGCTTGCCAAACCTCAGCCCCAGAA
TAAGGAAGCGGGAGTGAGGCAGAATCCTCAAGTGACAGCCAGGCCAGGAAGGAAGAGGAAGAGGAAC
AGGAAAACGGATGAGAAGGCTGGCAAGGGGCTCCTCTGGGAAGAGGAGGAAGCCGGGCTAGGGCAA
AGGGGGCTGGGGCCAGCAACCAAGAAGAGGATACAGTCATCACAGAGTGCAGAGTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	BC022910
Insert Size:	1179 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>BC022910, AAH22910</u>
RefSeq Size:	1866 bp
RefSeq ORF:	1178 bp
Locus ID:	18114
Cytogenetics:	10 39.72 cM
Gene Summary:	Plays a critical role in the generation of 28S rRNA.[UniProtKB/Swiss-Prot Function]