

Product datasheet for **MC218329**

Mtch1 (BC023135) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Mtch1 (BC023135) Mouse Untagged Clone
Tag: Tag Free
Symbol: Mtch1
Synonyms: C77849, PSAP
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >BC023135
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGGGAGCCGGAGCTGGTGCAGGAGCTCGGGCGGGCGCCGGCCGGAGTCGAGGCCCGCGCTCGGG
ACCCACCACCCGCGCACCAGCGCGCACCCCTCGCCATCCTCGGCCCGCGGCTCAGCCGTCGGCGCGCAGGAT
GGACGGCGGGCCCGGGCGCCCGGGCTCCGGGGACAACGCCCGACACCGAGGCGCTGTTCTGTGGCGCTG
GGCGCGGGCGGTGACGGCTCTCAGTCACCCGCTGCTCTACGTGAAGCTGCTGATCCAGGTGGGTGATGAGC
CGATGCCCCACCCCTTGGGACCAATGTGCTGGGAGGAAGTCCCTCTACCTGCCGAGCTTCTTACCTA
TGCCAAGTACATTGTGCAGGTGGATGGGAAGATAGGGCTCTTCCGGGGCCTGAGCCCCCGCCTTATGTCC
AACGCCTGTCCACTGTGACCCGCGGCAGCATGAAGAAGTTCCTCCAGATGAGATGGAGCAGGTTT
CCAACAAGGACGACATGAAGACCTCACTCAAGAAAGTTGTGAAGGAGACATCGTATGAGATGATGATGCA
GTGTGTATCGCGAATGCTGGCCATCCCTTACACGTGATCTCGATGCGATGCATGGTGCAGTTTGTGGGA
CGGGAGGCCAAGTACAGTGGTGTGCTGAGTTCTATTGGGAAGATCTTCAAGGAAGAGGGGCTGCTGGGAT
TCTTCGTTGGCTTAATCCCTCACCTCCTGGGCGATGTGGTTTTCTGTGGGCTGTAACCTGCTGGCCCA
CTTCATCAATGCCTACTTGGTGGACGACAGCGTGAGTGACACCCAGGGGGGCTGGGAAACGACCAAGAT
CCAGGTTCCAGTTTAGCCAGGCCCTGGCCATCCGGAGCTACACCAAGTTTGTGATGGGGATTGCAGTGA
GCATGCTGACCTACCCCTTCTGCTCGTTGGAGATCTCATGGCAGTGAACAAGTGTGGGCTGCGGGCTGG
ACTCCCTCCGATTTCCCTGTGTTCAAGTCTGGATCCACTGCTGGAAGTACCTGAGTGTGACAGGGCCAG
CTCTCCGCGGCTCCAGCCTGCTTTCCGCCGGGTGTCATCGGGGTCATGCTTTGCCCTGGAGTAA

ACGCGTACGCGGGCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
ACCN: BC023135



Insert Size:	1116 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC023135</u> , <u>AAH23135</u>
RefSeq Size:	1776 bp
RefSeq ORF:	1115 bp
Locus ID:	56462
Cytogenetics:	17 A3.3
Gene Summary:	Potential mitochondrial transporter. May play a role in apoptosis (By similarity). [UniProtKB/Swiss-Prot Function]