

Product datasheet for **MC218482**

Arhgap15 (BC034881) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Arhgap15 (BC034881) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Arhgap15
Synonyms:	5830480G12Rik
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

Fully Sequenced ORF: >BC034881
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAGAATCAAAAATGCCAACAGCCACCAAGACAGGCAGAGCCAAACTAAATCCATGATCCTCACCGATG
 CTGGGAAGGTCACTGAACCTATTTCCCGGCACAGAAGGAATCATTACAGCATGTGTTGAAAGATGTCAT
 TCCACCACTAGAACATCCAATGGTTGAAAAAGAAGGGTATCTTCAAAAAGCCAAAATTGCAGATGGAGGA
 AAGAACTAAGGAAAACTGGTCTACTTCTGGATTGTTCTTTCTGGTCGAAAAATTGAATTTTACAAGG
 ACTCCAAGCAGCAAGCTTCTTAACATGAAAACAGGCACAACGTGGAAAGTGTGGATTGTGTGGTGC
 ACATATAGAATGGGCCAAAGAAAAATCAAGCAGAAAGAGTGTCTTTCAGATCACAAACAGTGTCCAGAAAT
 GAGTTCCTTCTACAGTCAGATATTGACTTCCTCATATTGGATTGGTTCGAAGCTATCAAAAATGCAATTG
 ACAGATTGCCAAAGAATCCAAGTTGTGGTCCCTGGAGTTGTTCAATTTGCAGAGATCCTCAAGTTCTGA
 ACTGCCGAGTCACTGCCACATCGATAGAAAAGAACAGAAACAGAACACAGGAAGTCTTCATGTTTCGA
 CTGCACCACAGTCTTCTGATACAAGTGACAAGAATCGCGTGAAGAGCAGACTGAAGAAGTTCATCTCCA
 GAAGACCTTCTCTGAAAACCTTTCAGGAAAAGGGACTCATTAAAGATCAAATTTTGGCTCTCATCTGCA
 CACAGTGTGTGAACGAGAACATTCCACAGTTCCATGGTTTGTAAAGCAATGCATTGAAGCTGTTGAAAA
 AGAGGCCCTAGACGTTGATGGAATTTATCGAGTTAGTGGCAATCTTGCAACAATACAGAAGTTAAGATTTA
 TTGTCAACCAAGAAGAGAAGCTGAATTTGGATGACAGCCAGTGGGAGGACATCCACGTCTCACCGGAGC
 ACTGAAGATGTTTTCCGGGAGCTGTCTGAACCGCTCTCCCTTACAGTTCTTTGAGCGGTTTGTGGAG
 GCGATCAAAAAACAAGACAGCAATGAAAAAATTGAACTATGAGGTCTCTGGTAAAACGTCTCCCTCCAC
 CAAATCATGATACCATGAAAATCCTCTCAGACATCTAACCAAGATAGTGCCAAAGCCTCCAGAAATCT
 CATGTCCACCCAAAGCTTGGGGATTGTGTTTGGACCCACCTTCTGCGAGCGGAAAAATGAGTCAGGGAAT
 GTAGCGGTCCACATGGTATACAAAACAGATAGCAGAGTTCATGCTGACTGAGTACGATAAGATCTTCA
 GCTCAGAGGAAGACTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: BC034881

Insert Size: 1347 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC034881](#), [AAH34881](#)

RefSeq Size: 1724 bp

RefSeq ORF: 1346 bp

Locus ID: 76117

Cytogenetics: 2 B

Gene Summary: The protein encoded by this gene is a RAC GTPase-activating protein that is regulated through its PH domain and by recruitment to the membrane. The protein accelerates hydrolysis of guanosine triphosphate to guanosine diphosphate to repress Rac activity. Knock-out of Arhgap15 function demonstrates that this gene is required to regulate multiple functions in macrophages and neutrophils. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Sep 2014]