

Product datasheet for **MC218505**

Ntng1 (NM_030699) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ntng1 (NM_030699) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Ntng1
Synonyms:	A930010C08Rik; AI853992; Lmnt1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

Fully Sequenced ORF: >MC218505 representing NM_030699
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGATTTTGTCAAGATTCCTGTGTCGATCCATGCCCTGTGGGTGACAGTGTCTCTGTGATGCAGCCCTACC
 TTTTCGTGTGGGACATTATGATGTATGTAAAGCCTGATTTACACAGAAGAAGCAAAGTTGGGATTA
 CACAGCCTGCCAGCCGGAAATCCACGGACATGACCAAGTATCTGAAAGTAAAAGTGGACCCCTCCGGATATT
 ACCTGTGGAGACCCTCCAGAGTCTTCTGTGCAATGGGCAACCCTTACATGTGCAATAATGAGTGTGATG
 CGAGTACCCTGAACTGGCACACCCTCCTGAGCTGATGTTTGAAGGAAGACATCCCTCCACATT
 TTGGCAGTCTGCTACTTGAAGGAGTACCCAAACCTCTCCAGTTAACATCACTCTGTCTTGGAGCAAA
 ACCATTGAACTCACAGACAACATAGTTATTACCTTTGAATCGGGCGTCCAGACCAAATGATCCTAGAGA
 AATCTCTCGACTACGGACGAACATGGCAGCCCTATCAGTATTATGCCACAGACTGCCTCCATGCATTCCA
 CATGGACCCGAAATCCGTGAAGGATTTATCTCAGCACACGGTCTTGGAAATCATTTCGACGGAAGAGTAC
 TCCACTGGGTACTCCACGAATAGCAAAATATCCACTTCGAGATCAAAGCAGGTTTGCCTTTTCGCTG
 GACCTCGGCTACGAAATATGGCTTCCCTCTATGGACAGCTGGATACAACCAAGAAACTCAGAGATTTCTT
 CACTGTACAGACCTGAGGATCAGGCTGTTGAGACCCGCCGTTGGGGAAATATTTGTAGATGAACTACAT
 TTGGCAGGTTACTTTTATGCGATCTCAGACATAAAGGTGCGAGGAAGGTGCAAGTGAACCTGCATGCCA
 CTTTCGTGTTTGTATGACAACAGCAAACCTGACATGTGAATGTGAGCACAACACTACAGGTCCCGACTGTGG
 GAAATGCAAGAAGAACTACCAGGGCCGACCTTGGAGCCCCGGCTCATACCTCCCCATCCCCAAAGGCACC
 GCAAACACCTGTATCCCCAGCATTTCAGTATCGGTAAGTGTGAATGCTTCGGCCACTCCAATCGGTGCA
 GTTATATCGATCTGCTAAACACAGTCAATTCGCTGAGCTGTAACACAACACTAGAGGGCAGCAGTGTGA
 GTTATGACAGGCTGGGCTACTTCAGAAATGCTTCTGCACAGCTGGACGATGAGAATGTGTGCATAGAGTGT
 TATTGTAACCCCTTGGGCTCAATCCATGATCGTTGTAATGGCTCAGGATTTTGTGAGTGTAAAGACTGGAA
 CAACAGGGCCTAAATGTGATGAGTGTCTGCCAGGAAATCCTGGTACTACGGCTGTCAACCTAATGTCTG
 CGACAATGAGCTCCTGCACTGCCAGAATGGAGGGACCTGCCAGAACAATGTGCGCTGCGCGTGCCAGAC
 GCCTACACCGGCATCCTCTGTGAGAAGCTACGGTGCAGAGAGGGCGGCAGCTGTGGCTCCGAATCCGGCC
 AGGGAGACCCCCGCGGGCTCCCCAGCACTGCTGCTGCTGACCATGCTGCTGGGACTGCCGTTCCCT
 GGTGTTCTAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_030699

Insert Size: 1620 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_030699.2](#), [NP_109624.1](#)

RefSeq Size: 5347 bp

RefSeq ORF: 1620 bp

Locus ID: 80883

UniProt ID: [Q8R4G0](#)

Cytogenetics: 3 F3

Gene Summary: Involved in controlling patterning and neuronal circuit formation at the laminar, cellular, subcellular and synaptic levels. Promotes neurite outgrowth of both axons and dendrites. [UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (a) represents the longest transcript and encodes the longest isoform (a).