

## Product datasheet for **MC218526**

### Lin7a (BC105965) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Lin7a (BC105965) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Lin7a  
**Synonyms:** Veli1, LIN-7A, MALS-1, TIP-33, Veli  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >BC105965  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGCATGAAACGATTACTGTTAATGGCTGCCCTGAATTCGTCGAGGGCCACAGCAAAGGCAACAGTTG  
CAGCTTTTGCAGCCAGCGAAGGCCACTCCACCCTCGGGTAGTTGAACTGCCAAAGACTGATGAAGGCCT  
GGGTTTTAATGTGATGGGAGGGAAGGAACAGAATTCCTCGATTTACATCTCCCGCATCATCCTGGAGGG  
GTGGCTGAAAGACATGGAGGCCTCAAAGAGGGGACCAGCTGCTATCAGTGAACGGAGTGAGTGTGGAAG  
GGGAGCACCATGAGAAAGCTGTGGAACCTCAAGGCTGCAAAGGCAACAGTTGCAGCTTTTGCAGCCAG  
CGAAGGCCACTCCACCCTCGGGTAGTTGAACTGCCAAAGACTGATGAAGGCCTGGGTTTTAATGTGATG  
GGAGGGAAGGAACAGAATTCCTCGATTTACATCTCCCGCATCATCCTGGAGGGGTGGCTGAAAGACATG  
GAGGCCTCAAAGAGGGGACCAGCTGCTATCAGATACACCCCAAAGTCTGGAAGAGATGGAGGCTCGT  
TTCGAGAAGCTGCGGACAGCTCGGCGTCGACAGCAGCAGCTTGCTCATTGAGCAGCAGCAACAGCAGC  
AGCAACAACAACCACAACAAACCACATGTCATAGGCGCTCGAAGGAAAGCTAG

**ACGCGT**ACGCGGCGCTCGAGCAGAAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** BC105965  
**Insert Size:** 1029 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">BC105965</a>
<b>RefSeq Size:</b>	5500 bp
<b>RefSeq ORF:</b>	683 bp
<b>Locus ID:</b>	108030
<b>Cytogenetics:</b>	10 D1
<b>Gene Summary:</b>	Plays a role in establishing and maintaining the asymmetric distribution of channels and receptors at the plasma membrane of polarized cells. Forms membrane-associated multiprotein complexes that may regulate delivery and recycling of proteins to the correct membrane domains. The tripartite complex composed of LIN7 (LIN7A, LIN7B or LIN7C), CASK and APBA1 may have the potential to couple synaptic vesicle exocytosis to cell adhesion in brain. Ensures the proper localization of GRIN2B (subunit 2B of the NMDA receptor) to neuronal postsynaptic density and may function in localizing synaptic vesicles at synapses where it is recruited by beta-catenin and cadherin. Required to localize Kir2 channels, GABA transporter (SLC6A12) and EGFR/ERBB1, ERBB2, ERBB3 and ERBB4 to the basolateral membrane of epithelial cells.[UniProtKB/Swiss-Prot Function]