

Product datasheet for **MC218668**

Tcrb (BC030075) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids

Tag: Tag Free

Symbol: Tcrb

Synonyms: BV13S1

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC030075
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGCATCCAGACCCTCTGTTGTGTGATCTTTATGTTCTGATAGCAAATCACACAGATGCTGGAGTTA
CCCAGACACCCAGACATGAGGTGGCAGAGAAAGGACAAACAATAATCCTGAAGTGTGAGCCAGTTTCAGG
CCACAATGACCTTTTCTGGTACAGACAGACCAAGATACAGGGACTAGAGTTGCTGAGCTACTTCCGCAGC
AAGTCTCTTATGGAAGATGGTGGGGCTTTCAAGGATCGATTCAAAGCTGAGATGCTAAATTCATCCTTCT
CCTCTCTGAAGATTCAACCTACAGAACCCAAGGACTCAGCTGTGTATCTGTGTGCCAGCAGTTGGGGGCA
AGACACCCAGTACTTTGGGCCAGGCACTCGGCTCCTCGTGTAGAGGATCTGAGAAATGTGACTCCACCC
AAGGTCTCCTTGTGAGCCATCAAAGCAGAGATTGCAAACAAACAAAGGCTACCCTCGTGTGCTTGG
CCAGGGGCTTCTTCCCTGACCACGTGGAGCTGAGCTGGTGGTGAATGGCAAGGAGGTCCACAGTGGGT
CAGCACGGACCCTCAGGCCTACAAGGAGAGCAATTATAGCTACTGCCTGAGCAGCCGCCTGAGGGTCTCT
GCTACCTTCTGGCACAATCCTCGAAACCACTTCCGCTGCCAAGTGCAGTTCCATGGGCTTTAGAGGAGG
ACAAGTGGCCAGAGGGCTCACCAACCTGTCACACAGAATCAGTGCAGAGGCTGGGGCCGAGCAGA
CTGTGGAATCACTTCAGCATCTATCATCAGGGGGTCTGTCTGCAACCATCTCTATGAGATCCTACTG
GGGAAGGCCACCCATATGCTGTGCTGGTCACTGGCCTGGTGTGATGGCCATGGTCAAGAAAAAATT
CCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: BC030075



Insert Size:	447 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC030075
RefSeq Size:	1125 bp
RefSeq ORF:	914 bp
Locus ID:	269846
Cytogenetics:	6 19.71 cM 6 B1