

Product datasheet for **MC218714**

OTTMUSG00000016571 (BC098094) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: OTTMUSG00000016571 (BC098094) Mouse Untagged Clone
Tag: Tag Free
Symbol: OTTMUSG00000016571
Synonyms: OTTMUSG00000016571
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >BC098094
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGACTGTAACCAATGTGGTAAGGACTTTAGAACAAGGAATGACCTCAAATACATAAGCGAACACATAC
AGGAGAGAAAACCTATGAATGTAACCAATGTGGTAAAGCCTTTGCAAGCAGTAGTGACCTCGAAAAACAT
AAACGAACACATACAGGAGAGAAAACCTATGAATGTAACCAATGTGGTAAAGCCTTTGCAGAAAAGCAGTA
CTCTCAAATCCATAACCGAACACATACAGGAGAGAAAACCTATGAATGTAACCAATGTGGAAAAGCCTT
TGCAAGAAGGTGTGACCTCCAAAAACATAAACGAACTCATACAGGAGAGAAAACCTATGAATGTAACCCAC
TGTGGTAAAGCCTTTGCACAAAGCAGTTCTCTCCGAATCCATAAGCGAACACATACAGGAGAGAAAACCT
ATGAAAGTAAACAATGTGGTAAAGCCTTTGCACAAAGCAGTCATCTCCGAAAACATAACTGAACACATAC
AGGAGAGAGACAATGAATGTAAGTATGGTAAAGCCTTTGCAGGAAGCAGTGGTCTCCAATGCCAT
AAAAGATCACATACAGAAAGACACCATATGAATGTAACCAAGGTGGTAAACCTTTGAAGGAAGCAGTG
TTCTTGAATATCATAATCGAACACACAGCTGAAAAACCTATGAATTTAACCAATGTAGTAAAGCCTT
TGCAGGAAGGAGTGACCTTTAAAGAGAAAGGACCATACAGGAGAGAAAACCTATGTAGGTATCTGATGTT
ATAA

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTATCTCTGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: BC098094
Insert Size: 774 bp



[View online »](#)

| | |
|-------------------------------|---|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>BC098094</u> , <u>AAH98094</u> |
| RefSeq Size: | 3733 bp |
| RefSeq ORF: | 773 bp |
| Locus ID: | 433520 |
| Cytogenetics: | 2 |