

Product datasheet for **MC219006**

Tbc1d24 (NM_001163850) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tbc1d24 (NM_001163850) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Tbc1d24
Synonyms:	9630033P11; C530046L02Rik; mKIAA1171
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

Fully Sequenced ORF: >MC219006 representing NM_001163850
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGACCCCCAGGGTACAATTGCTTTGTGGATAAGGACAAGATGGATGCATCCATCCAGGACCTGGGGC
 CAAAGGAGCTGAACTGCACGGAGCTACAGGAGCTGAAGCAACTGGCACGGCAGGGTTACTGGGCTCAGAG
 CCACACCTGCGCGGAAAGTGTACCAGCGCTGATCCGGGACATCCCCTGCCGCACAGTACACCTGAT
 GCCAGCGTGTACAGTGACATTGTGGGAAGATTGTGGCAAGCACAGCAGCAGTAGTCTGCCCTTGCCCTG
 AGTTTGTAGACAACACTCAGGTGCCACCTACTGCCTGAACACACGGGGTGAAGGGGCTGTGCGCAAGAT
 CCTCTGTGATTGCCAACAGTCCCTGACATCTCCTTCTGCCCTGCCCTGCCTGCTGTGGTGGCCTTG
 CACTGCACTACAGCATCGATGAAGCTGAGTGTTTCGAAAAGCCTGCCGCATCTTATCCTGCAATGACC
 CCACCAAGAAGCTCATTGACCAGAGCTTCTGGCCTTTGAGTCTTCTGTATGACATTTGGGGACCTGGT
 GAACAAGTACTGCCAGGCAGCCATAAACTGATGGTGGCCGTGTCAGAGGACGTCTGCAGTCTACTCT
 GACTGGCAGCGATGGCTCTTTGGGGAGCTGCCCTCAATTAATTTGCCCCGTGTTTTGATGTCTTCCCTG
 TAGAAGGCTACAAAGTATTGTACCGAGTTGCTCTGGCCATCCTCAAGTCTTCCACAAGGTGAGAGCAGG
 ACAGCCCCCTTGAGTCAGACAATGTAAGCAAGACATCCGAATGTTTGTCAAGGACATTGCCAAGACAGTG
 TCCCCTGAGAACTACTAGAGAAGGCCTTTGCCATCCGCCTATTTCCCGAAAGGAGATCCAGCTTTGC
 AAATGGCCAATGAGAAAGCACTGAGGCAGAAGGGTATAACTGTCAAACAGAAGAGGCAGTTTGTGCACTT
 AGCTGTCCATGCAGAGAAGTCCACTCAGAGATTGTCAGCGTGAAGGAGATGAGAGACATTTGGTCTGG
 ATCCCTGAGCGCTTTGCCCTCTGCCAACCCCTCCTGCTTCTCATCACTGCAGCATGGTACAGCCTAA
 GCAGGTTCTATTTCCAGTGTGAAGGACATGAGCCACCCTCCTGCTCATCAAGACCACTAAAAGGAGGT
 CTGTGGGCTTACCTGTCAACAGACTGGAGCGAGAGGACTAAATTCGGAGGCAAGCTGGGCTTCTTTGGG
 ACTGGAGAATGCTTTGTTTTAGGCTGCAGCGGAAGTGCAGCGCTATGAGTGGGTGGTATCAAACACC
 CAGAGCTGACCAAGGCAACATCCCTCAAGTCTTCAGAGGCTGCAGGCAGTTCTTCCCTCATCAGCCACTG
 CTCTTCAGACCCCGCTGATCGGCTCTCCCTTTCTGGCTGCCCGGCACTTCAACCTGCCCTCCAAGACT
 GAGTCCATGTTTATGGCTGGAGGCAATGACTGCCTCATCATTGGGGGAGGAGGTGGTCAAGCACTTACG
 TTGATGGGACCTGAATCGAGGGCGCACTGGACACTGTGACACATTCAACAACCAGCCCTCTGCTCTGA
 GAACTTCTCATCGCAGCTGTGGAGGCTTGGGGCTTCCAAGACCCTGACACCGAA**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001163850

Insert Size: 1668 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001163850.1](#), [NP_001157322.1](#)

RefSeq Size: 7710 bp

RefSeq ORF: 1668 bp

Locus ID: 224617

UniProt ID: [Q3UUG6](#)

Cytogenetics: 17 A3.3

Gene Summary: May act as a GTPase-activating protein for Rab family protein(s) (PubMed:20727515). Involved in neuronal projections development, probably through a negative modulation of ARF6 function (PubMed:20727515).[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (5) lacks an alternate in-frame exon in the central coding region, compared to variant 1, resulting in an isoform (b) that is shorter than isoform a. Variants 4-8 encode the same isoform. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.