

Product datasheet for **MC219456**

Arid5a (NM_001172205) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Arid5a (NM_001172205) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Arid5a
Synonyms:	D430024K22Rik; Mrf1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF:

>MC219456 representing NM_001172205

Red=Cloning site Blue=ORF Orange=Stop codon

 TTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

 ATGGCAGCACCTCCGGCCAAAGGGAACACAGAGCAGTCAGAAGAAGGTGACCTCCCGCAGCTTCCTGTAT
 CCCCCAAGCCAGATGATGAGCAGAGCAGGAGCCAGAGCCCCACCCAGCTCCAGGACTCCCCTGAGGCAGG
 TGGGGAGCAGGAGGAGGAACAGGCCTTCCTGGTCAGCCTCTACAAGTTCATGAAGGAGCGACACAGCCCC
 ATCAGAGGGTGCCCCATCTTGGCTTCAAGCAGATTAACCTGTGGAAGATCTACAAGGCAGTGGAGAAGC
 TGGGGGCCTATGAGCTGGTGACAGGCCGCCCTCTGGAAGAAGTGTATGATGAACTTGGCGGTAGCCC
 AGGCAGCACCAAGTGGGCCACATGCACAGCGGCCACTATGAGAGGCTGGTCTCCCATATGTGCGGCAT
 CTGAAGGGGGAGGACGACAAGCCACTGCCTCTACCAAGCCCAGGAAGCAATAACAAGATGGCCAAGGAGC
 TGAGGGGAGACGATGGGACCACTGAGAAGCTGAAGAAGGCCAAGGACTCAGAGGAGAGGCGGGTGGAGCA
 GACCACGCCAGAAAGACCAATCAGATGCCACTGGCCAGACACAGCTTCCCTGCCAGGGATCCTCGAGG
 GACAGCACAGAACAGCTGGGCCAGTATCTGGACCCTCTCCACCACTCACGGGTGCTAGTAGCTGCCCTG
 AGGCCTACAAGCGGCTCTTGTCAGCTTTTACTGCAAAGGGGCGCATGGCATCATGTCAACCACTGGCCAA
 AAGAACTCTTGGCCCAGGTCAGCAAGGCAGAGGCCCTTGCAAGTCCCAAGAGGGCTGTGCCATGGA
 GCAAGGAGCCCCAACAGGACATTCAAGACAGTCCCCAGAACCTAAGAGGGCCGGCTGAGAACTCTGAAC
 ACCAGTAACCCCCCGGAAGGATTGCAGGCCCTGGTGGGAGCACCAGGATGGAGGCCAAGTGGGCC
 CTGCCCTACAGCCCCATGTTCTCAGGCTGTTTTATGCGTACCCACCGAGGTGCTGAACTGTGAGC
 CAGCACCTAGGACTTCTTCTCCGGCCTTAAAGACAGGGTGTGTTGGGACCACCTGGTAAAGAAGAAG
 GTCCGACAACCAAGAGTCCCCTCTGGTGTGGGTGGGATGCCAACCACCCCTCTGCATTCCATAAAGG
 CAGCACAAAGAAAAAGATTCTACCCCAACCCAAAGCCTGCTGGGTGTCTCCCATGGCCAAGGTCCCT
 ACTGAGAGGCCTGGAGCCCCATCCCTCATCCAGTAGCCAGGTCTTGCAAGTAAAGCGCGGCTTGAAG
 AAGAGGATTGCTCATGGTGGCAAGAACTGAGGGCAGTGTCTCCCTTCTGAAGGAGGTGGATTCCAA
 GGAGACTGGGGCAAGCCTGCAGCCCCTGGCTTGCTGTATCCTGTCTACTGGGCCAACCCGGGGCCC
 ACTCCTCCAGAGGCCTACAGGGGCACCATGCTGCGGTGTCTCTAACTTCAACGGTAGCGCAGACCCTC
 TGAAGGGCCAGGCCTCACTCCCCTTCAAGCCCCTGGTCATCCCTGCTTCCAGCCACCTTCTGGCTAC
 AACAGGCTCCTCACCTATGGCTGCCAGCCTGATGATTTCCCTCCACGCCCTATGACGCTGTCTACGC
 AACAGACTGGGTCCAGCTTCGTCTGCCTGGCACATGCCACCCGTACAACCTATGCGGCACCTCACTTCT
 TCCACCTCAACACCAAACT**GTAG**
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:

SgfI-MluI

ACCN:

NM_001172205

Insert Size:

1773 bp

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001172205.1 , NP_001165676.1
RefSeq Size:	5699 bp
RefSeq ORF:	1773 bp
Locus ID:	214855
UniProt ID:	Q3U108
Cytogenetics:	1 B

Gene Summary:

DNA-binding protein that may regulate transcription and act as a repressor by binding to AT-rich stretches in the promoter region of target genes (By similarity). May positively regulate chondrocyte-specific transcription such as of COL2A1 in collaboration with SOX9 and positively regulate histone H3 acetylation at chondrocyte-specific genes. May stimulate early-stage chondrocyte differentiation and inhibit later stage differentiation (PubMed:21346191). Can repress ESR1-mediated transcriptional activation; proposed to act as corepressor for selective nuclear hormone receptors (By similarity). As RNA-binding protein involved in the regulation of inflammatory response by stabilizing selective inflammation-related mRNAs, such as IL6, STAT3 and TBX21. Binds to stem loop structures located in the 3' UTRs of IL6, STAT3 and TBX21 mRNAs; at least for STAT3 prevents binding of ZC3H12A to the mRNA stem loop structure thus inhibiting its degradation activity. Contributes to elevated IL6 levels possibly implicated in autoimmunity processes. IL6-dependent stabilization of STAT3 mRNA may promote differentiation of naive CD4+ T-cells into T-helper Th17 cells (PubMed:23676272, PubMed:27022145). In CD4+ T-cells may also inhibit RORC-induced Th17 cell differentiation independently of IL6 signaling (PubMed:24782182). Stabilization of TBX21 mRNA contributes to elevated interferon-gamma secretion in Th1 cells possibly implicated in the establishment of septic shock (PubMed:27671645).[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (1) represents the predominant transcript and encodes isoform 1.