

Product datasheet for **MC219473**

Fbrsl1 (NM_028596) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fbrsl1 (NM_028596) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Fbrsl1
Synonyms:	2410025L10Rik; AI467568; Gm1049; Gm29766; mKIAA1545
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF:	>MC219473 representing NM_028596 Red=Cloning site Blue=ORF Orange=Stop codon TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC GCCCGGATCGC ATGTTGCTGCCCCCAACACTGCCCCGCCCCGGCGTTGCCCGCCAGCAGCCTGGTCATCCCAGGAC ACCCCGCGATCACGAACTGCTCAGGCAAGAGCTGAACACGCGGTTTCTGGTGCAGAACGCGGAGCGGCC TGGCGCTCCCTGGGCGCGGGGCTCTGCTGCGGGCCGAGTTCACCAGCACACGCCAGCACACAC CAGCACACACACAGCACCAACACACATTTGCCCCCTTGCAGGGCTGCCCCGACGCCGCTTATGCCGC CCACCGCACCTCGCCGTTTGACAAGTTCTACCCAAGCTGGACAGCCCTACTTCCGACATTCCAACAC TTCAAACCTCAATCGAAATAACAGGCCGGTCGAGTGCAGTTCACACTTTCCTGCAGAAAGCCCTGGGGTG TCCGATCCGTACCGACAGCAGTCAGGATGCAGCTGGACCCCCACAAGCTAGAGGTTGGCAGAAAGCTGG ACCTGTTACGCGGCCCCCGCCCCAGGCATGTTCTCTGGGTTCCACTACCCACAGGACCTGGCCCCGCC CCTCTTCTCCAGCACAGGTGCCACCCATCCAAGTCCAAACCCATTCGACCCCTCAGCTCATCATGGCAGC TTCCTGCCCACTGGTCACCTGACAGACCCTTTCAGCAGATCCAGCACTTTCGGAGGCCTGGGGAGCCTGA GCAACCACGCGTTTGGAGGTCTGGGCAGCCATGCACTGACACCAGGAAGCAGCATCTTGGCCCCAAGGA GGGCACAGCTTTGCTCGGCCTGCCAGCCCCCATGATGCCTGGAACCGGCTGCACCGGTGCCACCTTCC TTTCTACCCACCGCCATGGCCCAAGCCAGTGGACACAGAGCGCCTCTCAGCCCTGACCAACCATGACC GAGAGCCTGACAAGGGCCAGGAAGAGCATGAGAGGGACCTCTAGAGAAGACGCGCCTGCTGAGCCACGG GTCGCCTGTAGCACCAGTGGGCCACCCATCAGCAGCCTCCTACAGCGGGGCTCCGGTGAGCTGAGCCGC CCTGGGGTCTCACCAGCGTGAGGCTGAGACTCGGGTCAAGGAGAGCCACTCACCTAGCAAAGAGGATG GCGCCAAACTGGCTGCGCGACCTCCGTACCCCTACTGCAAGGCAGCCCTGGGCGACTGCCTACGCCTGGC AGGCCTACTGGGTCGGGAGTCCGGGAAGCAGACGCTGATCGGCCCCAGGGCGACGTGAAGGTCAAGGAG GAACGCAGGGAAGACGAGACACCGGACCTCACGGGCGTCAGCGTGCCTCCACAGCTACACATATGGGCC TGGGCACAGCAGGCTTCGTGTGGGATCCGCCGCGGACACCTACCGGGGCTGGAATACCGCGCGGCTC CTTCCAGCCCCGCGCCCGCTGCGGGACCTTCTCCCTTCGAGCTCTCGGAGCGCGCTACCGGGACCGC GAGCCTCAGACTACAGCCGGAGCGCCTGCAGGAGTGCGCCGGGAGGAGCTGGAGCGCGCGCGCTGC CCGATTCGACGCCGCGGTTGCTGCCGGCGTGCACTACCCACGCCTGGCTCCTGGGGCGCTGCACGG GGCGCTCCTGGCAGCACCCCGCGGCCGCGCGGTGCTCGGTGCGCCTCCTCCACTGGTGAAGTGCAGGGC GGGGTCCCCACGCCCGCGCATCCCAGGACGACGCTGGGGCCTGGCGAGGCACGGGAGTACTCGCCCC CGTGCAACCCTCCCGAGGTGGAGGCGCGTAG AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC TGGATTACAAGGATGACGACGATAAGGTTTAA
Restriction Sites:	Sgfl-RsrII
ACCN:	NM_028596
Insert Size:	1782 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_028596.2, NP_082872.2

RefSeq Size: 3781 bp

RefSeq ORF: 1782 bp

Locus ID: 381668

Cytogenetics: 5 F