

## Product datasheet for **MC219511**

### **Nr4a2 (NM\_013613) Mouse Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Nr4a2 (NM_013613) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Nr4a2
Synonyms:	HZF-3; NOT; Nurr1; RNR-1; TINOR; TINUR
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_013613, the custom clone sequence may differ by one or more nucleotides

```

ATGCCTTGTTGTTTCAGGCGCAGTATGGTCTCGCCTCAAGGAGCCAGCCCCGCTTCTCAGAGCTACAGTT
ACCACTCTTCGGGAGAATACAGCTCCGATTTCTAACTCCAGAGTTTGCAAGTTTAGCATGGACCTCAC
CAACACTGAAATTACTGCCACCCTTCTCTCCAGCTTACAGTACCTTTATGGACAACACTACAGCACAGGC
TACGACGTCAAGCCACCTTGCTTGACCAAATGCCCTGTCCGGACAGCAGTCTCCATTAAGGTAGAAG
ACATTCAGATGCACAACCTACCAGCAACACAGCCACCTGCCCTCAGTCCGAGGAGATGATGCCACACAG
CGGGTCGGTTTACTACAAGCCCTCTTCGCCCCGACACCCAGCACCCGAGCTTCCAGGTGCAGCATAGC
CCGATGTGGGACGATCCGGGCTCCCTTCAAACTTCCACCAGAACTACGTGGCCACTACGCATATGATCG
AGCAGAGGAAGACACCTGTCTCCCGCTGTCACTTCTCTTTAAGCAGTCCGCCCCGGGCACTCTGT
GTCTAGCTGCCAGATGCGCTTCGACGGGCTCTGCACGTCCCATGAACCCGAGCCCGGGCAGCCAC
CACGTAGTGGATGGCAGACCTTCGCCGTGCCAACCCTTTCGAAGCCGGCATCCATGGGCTTCCCGG
GCCTGCAGATCGGCCACGCATCGCAGTTGCTTGACACGCAGGTGCCCTCGCCGCGTCCCGGGCTCTCC
CTCCAATGAGGGTCTGTGCGCTGTTGCGGTGACAACCGGCCTGTGAGCACTACGGTGTTCGCACTTGT
GAGGGCTGCAAAGGTTTCTTTAAGCGCACGGTGCAAAAAACCGCAAATATGTGTGTTTAGCAAATAAAA
ACTGCCAGTGGACAAGCGCCGCGAAACCGTTGTGAGTACTGTGCGTTTTCAGAAAGTGCCTAGCTGTTGG
GATGGTTAAAGAAGTGGTTCGCACGGACAGTTTAAAGGCGGAGAGGTCGTTTACCCTCGAAGCCGAAG
AGCCACAGGATCCCTCTCCCCCTCACCTCCGGTGAAGTGTGATCAGTCCCTCGTCAGAGCCACGTCG
ATTCCAATCCGGCAATGACCAGCCTGGACTATCCAGGTTCCAGGCAAACCTGACTATCAGATGAGTGG
AGATGATACCCAAACATATCCAGCAGTTCTACGATCTCTGACCGGCTCTATGGAGATCATCAGAGGGTGG
GCAGAGAAGATCCCTGGCTTTGCTGACCTGCCAAAGCCGACAGGACCTGCTTTTGAATCAGCTTTCT
TAGAATTTATTTGTTCTGCGCTTAGCATACAGGTCCAACCCAGTGGAGGTAACCTATCTTTTGAATGG
GGTGGTCTTGACAGGTTGCAATGCGTGCCTGGCTTTGGGAATGGATTGATTCCATTGTTGAATTTCC
TCCAATTCAGAAATGAAATCGACATTTCTGCCTTCTCTGCATTGCTGCCCTGGCTATGGTACACAG
AGAGACACGGGCTCAAGGAACCAAGAGAGTGGAAAGCTACAAAACAAAATTGTAATTTGCTTAAAGA
CCATGTGACTTTCAATAATGGGGGTTTGAACCGACCAACTACCTGTCTAACTGTTGGGAAGCTGCCA
GAACTCCGCACCCTTTGCACACAGGGCCTCCAGCGCATTTTCTACCTGAAATTGGAAGACTTGTACCAC
CACCAGCAATAATTGACAACTTTTCTGGACACCTTACCTTTCTAA
    
```

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_013613

**Insert Size:** 1797 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">BC137715</a> , <a href="#">AAI37716</a>
<b>RefSeq Size:</b>	2642 bp
<b>RefSeq ORF:</b>	1797 bp
<b>Locus ID:</b>	18227
<b>UniProt ID:</b>	<a href="#">Q06219</a>
<b>Cytogenetics:</b>	2 31.66 cM
<b>Gene Summary:</b>	<p>Transcriptional regulator which is important for the differentiation and maintenance of meso-diencephalic dopaminergic (mdDA) neurons during development. It is crucial for expression of a set of genes such as SLC6A3, SLC18A2, TH and DRD2 which are essential for development of mdDA neurons.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longer transcript.</p>