

Product datasheet for **MC219751**

Lpp (NM_001145952) Mouse Untagged Clone

Product data:

| | |
|--------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | Lpp (NM_001145952) Mouse Untagged Clone |
| Tag: | Tag Free |
| Symbol: | Lpp |
| Synonyms: | 9430020K16Rik; AA959454; AU024130; B130055L10Rik; C79715; D630048H16 |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Cell Selection: | Neomycin |



[View online »](#)

Fully Sequenced ORF: >MC219751 representing NM_001145952
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTCTCACCCATCTTGGCTGCCACCCAAAAGCACTGGTGAAGCCTCTTGGCCACGTGCCTGCACGGATGG
 AGACCACCCATTCTTTGGGAACCCAGTATCTCAGTGTCAACACAACAGCCACCCAAGAAGTATGCTCC
 AGTTGTTGCTCCGAAACCAAGTATAACCCATACAAGCAGCCTGGAGGAGAGGGTATCTCTCCACCCG
 CCACCACCACTTAGAGGATCCTGGTACCATCCACCTGGTCCCGGACACTTCCCTCCTCCTCCACCAC
 TGGATGAAGGGGCTTTCAAAGTGCAGCAGGAAATCCTGGAGGCAAGACCCTGGAGGAAAGGCGCTAG
 TCTAGATGCAGAGATCGACTCCTTGACCAGCATCTTGGCTGATTTGGAGTGCAGCTCACCTACAAGCCA
 CGTCCCCACAGGGCTCTGCTAGTTCAATAGCCTCCTCCAGTTTCTACTCCTGTCACAGGACATAAGA
 GAATGGTCATTCCACAGCAGCCACCCCTGACAGCAACCAAGAAGTCCGCAACAAAACCACAGCCTGCACC
 CCAGGCTGCACCCATCCCAGTACTCCAATTGGAACACTCAAACCCAGCCTCAGCCAGTCCAGCCTCC
 TATACTACGGCATCCACATCTTCAAGGCCACCTTTAATGTGCAAGTGAAGTCAAGTCCAGCCAGCCAC
 ATTACATGGCTGGTCCATCCTCTGGACAAATATATGGTCCAGGGCCTCGTGGCTATAACAACCAGCCAGT
 TCCTGTATCGGGACAGTGTCTCCTCCTCCAACCTGTGTTGGTACCGACTATGCATATATCCACCATCT
 GGACATCCGCCTGAGTCTGGATATGGGTATACCTCAACCAAGGACGCTATTATGAACCTTACTATGCAG
 CAGGTCCTAGCTATGGGGCCGAAGTGAAGGTGATACTGCTTATGGGCAACAGGTTCAACCTAATACCTG
 GAAGCGGGAAGCTGCGTATGCCCTCCTGCAAGCGGGAATCAGAACCATCCTGGCATGTATCCAGTCAGT
 GGCCCCAAGAAGACCTATACAGACCCTGTTTCAGCGCCCTGTGCACCACCACTGCAGCCAAAGGGTG
 GATACCTGGGCCTATGGGGCCTCCTTCCATCCCCCATCATTCCGCCCAGAAGATGAGCTCGAACACTT
 GACCAAGAAGATGCTGTATGACATGGAGAACCACCCGCTGACGATTACTTTGGCCGTTGTGCCCGCTGT
 GGGGAAAATGTTGTTGGAGAAGGAACAGGATGACTGCCATGGATCAGGCTTCCACGTGGATTGCTTCA
 CCTGCATTGTCTGTGATGTCAAGCTCCGAGGGCAGCCTTTCTATGCCGTGGAGAAGAAGGCATACTGCGA
 GCCCTGTACATTAATACCCTGGAGCAGTGCAGTGTGTGTTGGAAGCCCATCATGGAGCGGATCCTCCGT
 GCCACCGGGAAGGCTACCACCCCACTGCTTACCTGCGTGTGTGCCACCGTAGCCTGGACGGCATCC
 CATTACTGTGGATGCTTGGGCTCATTCACTGCATTGAGGACTCCACAAGAAATTTGCCCGCGCTG
 TTCTGTGTGCAAGAACCTATCATGCCAGCCCCAGGCCAGGAGACTGTCCGATTGTGGCTCTGGAC
 CGTGATTTCCATGTGCACTGCTATCGCTGTGAGGACTGTGGTGGTCTTCTGTCCGAGGAGACAACCAAG
 GCTGCTACCCCTTGATGGACACATTCTCTGTAAGACCTGCAACTCGGCTCGCATCAGGGTCTTGACGGC
 CAAGGCAAGCACTGACCTG**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI

ACCN: NM_001145952

Insert Size: 1842 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001145952.1](#), [NP_001139424.1](#)

RefSeq Size: 15462 bp

RefSeq ORF: 1842 bp

Locus ID: 210126

UniProt ID: [Q8BFW7](#)

Cytogenetics: 16 B1

Gene Summary: May play a structural role at sites of cell adhesion in maintaining cell shape and motility. In addition to these structural functions, it may also be implicated in signaling events and activation of gene transcription. May be involved in signal transduction from cell adhesion sites to the nucleus allowing successful integration of signals arising from soluble factors and cell-cell adhesion sites. Also suggested to serve as a scaffold protein upon which distinct protein complexes are assembled in the cytoplasm and in the nucleus (By similarity).

[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Both variants 1 and 2 encode the same isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.