

Product datasheet for **MC219786**

Keap1 (NM_001110307) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Keap1 (NM_001110307) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Keap1
Synonyms:	INRF2; mKIAA0132
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

Fully Sequenced ORF: >MC219786 representing NM_001110307
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCAGCCGAACCCAAGCTTAGCGGGCTCCCGCAGCAGCCAGTTCCTGCCCTGTGGTCAAAGTGCC
 CCGAGGGGGCCGGGACGCAGTGATGTATGCCTCCACGGAGTGCAAGGCAGAGGTGACGCCCTCGCAGGA
 CGGTAACCGAACCTTCAGCTACACACTAGAGGATCACACCAAGCAGGCTTTTGGCGTCATGAACGAGCTT
 CGCCTGAGCCAGCAACTCTGTGACGTGACCCTGCAGGTCAAATATGAGGACATCCAGCTGCCAATTCA
 TGGCTCACAAAGTGGTGTGGCCTCTCCAGCCAGTCTTTAAAGCCATGTTACCAACGGGCTTCGGGA
 GCAGGGCATGGAGGTGGTGTCCATCGAAGGCATCCACCCTAAGGTCATGGAAAGGCTTATTGAGTTCGCC
 TACACGGCCTCCATCTCCGTGGCGGAGAAGTGTGCTGCACGTGATGAACGGGGCGGTCATGTACCAGA
 TTGACAGCGTGGTTTCAGCCTGCAGCGACTTCTCGTGCAGCAGCTGGACCCAGCAACGCCATTGGCAT
 CGCCAACTTCGCGGAGCAGATCGGCTGCACTGAACTGCACCAGCGTGCCCGGAGTATATCTACATGCAC
 TTCGGGGAGGTGGCAAGCAGGAGGAGTCTTCAACCTGTCACTGCCAGCTGGCCACGCTCATCAGCC
 GGGATGATCTGAACGTACGCTGCGAGTCCGAGGTGTTCCACGCGTGCATCGACTGGGTCAAATACGACTG
 CCCGCAGCGGCGCTTCTACGTGCAGGCACTGCTGCGGGCCGTGCGCTGCCATGCGCTCACGCCGCGCTT
 CTGCAGACGCAGCTGCAGAAAGTGTGAGATCCTGCAGGCCGACGCGCGTGAAGGACTACCTGGTGCAGA
 TATTCCAGGAGCTCACGTGCACAAGCCACGCAGGCAAGTGCCTGCCGCGCGCCAAAGTGGGCCGCT
 CATCTACACAGCGGGCGGTTACTTCCGACAGTCGCTCAGCTACCTGGAGGCTACAACCCGAGCAATGGC
 TCCTGGCTGCGCTGGCCGATCTACAGGTGCCGCGCAGTGGGCTGGCAGGCTGCGTGGTGGGTTGGCTGC
 TATACGCTGTGGGCGGCCGCAACAACCTCCGGATGGCAACACTGACTCCAGCGCCCTGGACTGCTACAA
 CCCCATGACCAACCAGTGGTCCGCTGTGCTCTATGAGCGTGCCACGCAACCGCATCGGGGTGGGGGTC
 ATAGATGGCCACATCTACGCAGTCGGGGTTCACGCGTGCATCCACCACAGCAGCGTGGAGAGATATG
 AGCCAGAGCGGGACGAGTGGCATCTAGTCGCGCAATGTTGACACGGAGGATTGGCGTGGGCGTGGCAGT
 GCTCAACCGCTTGTGTATGCAGTGGGGGCTTTGACGGGACTAACCGGCTTAACTCCGAGAAATGTTAC
 TATCCAGAGAGGAATGAGTGGCGGATGATCACACCGATGAATACCATCCGAGCGGGGCGGGGCTGCGG
 TGCTGCACAACGTATCTATGCAGCAGGGGCTACGATGGCAGGACCAGTTGAACAGTGTGGAGCGCTA
 CGAGTGGAGACAGAGACCTGGACTTTCGTAGCCCCATGAGGCATCACCGTAGTGCCTGGGGATTACT
 GTGCACCAGGGCAAGATCTACGTCTCGGAGGCTATGATGGCCCACTTTCTGGACAGTGTGGAATGCT
 ATGACCCGGACAGTGATACCTGGAGTGAGGTGACCCGATGACATCTGGCCGACGCGGGTGGGTGTGGC
 CGTCACCATGGAACCCTGTCGGAAGCAAATTGATCAACAAAACGTACCTGCT**GTA**

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

- Restriction Sites:** SgfI-RsrII
- ACCN:** NM_001110307
- Insert Size:** 1875 bp
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001110307.1](#), [NP_001103777.1](#)

RefSeq Size: 3151 bp

RefSeq ORF: 1875 bp

Locus ID: 50868

UniProt ID: [Q9Z2X8](#)

Cytogenetics: 9 A3

Gene Summary: Acts as a substrate adapter protein for the E3 ubiquitin ligase complex formed by CUL3 and RBX1 and targets NFE2L2/NRF2 for ubiquitination and degradation by the proteasome, thus resulting in the suppression of its transcriptional activity and the repression of antioxidant response element-mediated detoxifying enzyme gene expression. Retains NFE2L2/NRF2 and may also retain BPTF in the cytosol. Targets PGAM5 for ubiquitination and degradation by the proteasome (By similarity).[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (4) uses a different splice site in the 5' UTR, compared to variant 1. Variants 1, 2, 3, and 4 all encode the same protein.