

Product datasheet for **MC219905**

Pde1c (NM_001025568) Mouse Untagged Clone

Product data:

| | |
|--------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | Pde1c (NM_001025568) Mouse Untagged Clone |
| Tag: | Tag Free |
| Symbol: | Pde1c |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Cell Selection: | Neomycin |



[View online »](#)

Fully Sequenced ORF: >MC219905 representing NM_001025568
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGTCTCAACCAAGGAAATGAAGAATTCGAGAGCAACTCTCTGAAACACCTGCAACCTGAGCAGA
TCGAGAAAATCTGGCTTCGGCTCCGCGGGCTGAGGAAATATAAGAAAACGTCCCAGAGGTTACGGTCTTT
GGTCAAGCAATTAGAACGAGGAGAGGCCTCTGTGGTAGACCTTAAGAAAGATCTGGAATATGACGCCACA
GTGCTTGAATCTGTATATATTGATGAAACAAGCGACTGCTGGATACAGAGGATGAGCTTAGTGACATTC
AATCAGATGCTGTGCCTCAGAAGTTCGGGATTGGCTGGCCTCCACCTTCACACGCAGATGGGGATGAT
GCTTAGGAGAAGTGATGAGAAGCCAGGTTCAAGAGCATCGTCCATGCGGTGCAAGCTGGGATATTTGTG
GAAAGAATGTACAGACGGACATCAAACATGGTTGGGCTGAGCTATCCACCGGCTGTAATCGATGCATTGA
AGGATGTGGATACGTGGTCTTCGATGTCTTTCCCTCAATGAGGCCAGTGGAGATCATGCACTGAAGTT
CATTTTCTATGAATTACTCACACGTTATGACCTGATCAGCCGTTTTAAGATACCGATTTCTGCATTGTCT
TCATTTGTGGAGGCCCTCGAAGTGGGTTACAGCAAGCACAAAAATCCTTACCATAACCTGATGCATGCAG
CTGACGTCACACAGACTGTGCATTACCTCCTTTATAAGACAGGAGTAGCAAACCTGGCTGACAGAGCTGGA
GATCTTTGCAATAATCTTTTCGGCTGCCATCCATGACTATGAACATACTGGAACACAAAACATTTCCAC
ATTCAGACTCGGTGAGATCCAGCTATCCTGTACAATGACAGATCTGTCTGGAGAACCACCACCTGAGTG
CAGCCTACCGCCTTCTGCAGGAAGACGAAGAGATGAATATTCTGGTTAACCTCTCAAAGGATGACTGGAG
GGAGTTTCGAACCTTTGGTTATTGAGATGGTAATGGCCACAGATATGTCTGTCAATTTCCAGCAAATCAA
GCCATGAAGACAGCCCTGCAGCAACCAGAAGCAATTGAGAAGCCGAAAGCCTTATCCCTCATGTTACACA
CAGCAGACATCAGTCATCCTGCGAAAGCATGGGACCTGCACCAGCTGGACCATGTCTCTCCTGGAGGA
GTTCTTTAGACAGGGTGACAGAGAAGCAGAGCTGGGGCTGCCATTTTCTCCTTTGTGACAGAAAAGTCA
ACCATGGTTGCTCAGTCAACAAGTGGGTTTTATTGACTTATTGTGGAGCCACCTTCACTGTGCTCACGG
ACATGACCGAGAAGATTGTGAGTCCATTAATTGACGAAAGCTCCCAGACTGGTGGGACAGGGCAGAGGAG
ATCAAGTTTGAACAGCATCAACTCATCAGATGCAAAGCGATCCGGTGTCAAGAGTTCTGGGTGAGATGGA
AGTGCTCCCATCAACAATTCTGTCACTTCTGTTGACTATAAGAGTTTTAAAGCCACTTGGACTGAGGTGG
TGCAGATCAATCGGGAGCGGTGGCGAGCCAAGGTACCCAAAGAAGAAAAAGCCAAGAAGGAAAGCTGAAGA
GAAGGCTCGCCTGGCTGTGAGGAAAAGCAAAGGAAATGGAAGCCAAAAGCCAAGCTGAACAAGGCACA
ACCAGCAAAGGCGAGAAAAGACATCAGGAGAAGCCAAAAGTCAAGTCAATGGAACACGCAAGGGTGATA
ACCCTCGTGGGAAGAAGCTCAAAGGAGAAAAGGCAGGCGAAAAGCAGCAGAACGGTGACTTGAAAGACGG
TAAAAATAAGGCAGACAAGAAGGATCACTCCAACACCGGAAATGAGTCAAAGAAAACAGATGATCCAGAG
GAGTAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001025568

Insert Size: 1896 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

| | |
|-------------------------------|---|
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_001025568.2 , NP_001020739.1 |
| RefSeq Size: | 8864 bp |
| RefSeq ORF: | 1896 bp |
| Locus ID: | 18575 |
| UniProt ID: | Q64338 |
| Cytogenetics: | 6 27.65 cM |
| Gene Summary: | <p>Calmodulin-dependent cyclic nucleotide phosphodiesterase with a dual-specificity for cAMP and cGMP, which are key regulators of many important physiological processes (PubMed:8810348). Exhibits high affinity for both cAMP and cGMP (By similarity). Modulates the amplitude and duration of the cAMP signal in sensory cilia in response to odorant stimulation, hence contributing to the generation of action potentials. Regulates smooth muscle cell proliferation. Regulates the stability of growth factor receptors, including PDGFRB (Probable).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) represents the longest transcript. Variants 2 and 8 encode the same isoform (b). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p> |