

Product datasheet for **MC220222**

Limd1 (NM_013860) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Limd1 (NM_013860) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Limd1
Synonyms:	AW822033; D9Erttd192e
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >MC220222 representing NM_013860
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGATAAATACGATGATCTGGCCTTGAGGCCAGTAAGTTCATCGAGGACCTGAACATGTATGAGGCCT
 CCAAGGATGGGCTCTCCGGGTGGACAAGGGTGTGGCAACAACCCAGAATTTGAGGAACTCGAAGGGT
 GTTCGCGACCAAGATGGCCAAAATCCACCTCCAGCAGCAGCAACAGCAGCAGCTCTACAGGAGGAGGCC
 CTACCTAGGGCAGGCAGAAGCCCGGTCAACGGTGGGAACCGTCAGGGTGCAGAGCGCAAGCTGGCTGCAG
 ATGGGGCCGCTAAGCCTCCTCTGGCTGTGCCAACAGTGGCACCTGGATTAGCTACTACCACTGCGGCCGC
 GCAGCCCTCGTACCCATCTCAGGAGCAGAGGATCAGGCCATCTGCCCATGGTGAAGGCCTGGCAGTCAG
 AACTGTGGTTCAGGGAGGGGCTGTGAGTTCAGAGACCTGCTTTACACGGCCTGAGTCCCTCCTGTG
 AAGACCTTCCTGCCTCACTCACGGAGACTATTATGACAATTTCTTTTGGAAGCCACAGTGGGGTGA
 TAAACCTGAAGGGTGGCCAGCGTGAGTCTGGGTGTAGGAAGTGGATGGCCTGGTTGCCAGGGAATGAC
 TCCACATTGCCCAAATCCTGTGGAGACCATCACCTTACCAGCCACAGCTCTCCACAGTGTGCTCTGGCA
 GGTCTTTTGAAAGTGGCATCAGTGGCCAGGACGGTGGCATTGGTGGCCATAGCAGCGAGAAGCCAACAGG
 CCTTTGGTCCACTGCCTCCTCTCAGCGAGTGAACCTCGGTTTTCTTCCATGGGCTTGGAGAATGGGACC
 TCCGCTCAACCAAGGGCACAACGTGTTTCAGCACCGATGGTCCCTAGCAGCGCCAGCCAAGGAGCTTGTC
 CGAAAAGAGATTAGGTCTGGGATATGAGGCTTCAGGCAGGGTCTTCAAACCCCTTGTGGACACTCAGCC
 TTGGCTGCAGGATGGGCCAAGTCTTACCTCTCAGTCTTCTGCTCCATTATCCTCGACAGCTGGCAAGGAC
 AGTACCCAGCCAGGTATGACCACCGGGCTGGATCCTAAGTTTGGATGCGTGGAGTCTGGCACTAGTCCCA
 AGCCCAGCCCCACAGTAACGTCCATCCGTAATGTCCACTCCATCTGAGTTATCTTGAAGAGAGTTT
 TCCAGCTGGTCCACTGACAGTAGTCTGGAACCTGTGCTCCAGGGAGTCCCACCCCTCCAGGGTGAGA
 TTGCCCTGCCAGACCCTCGCACCGGGCCCTGAGCTTGGACCCTCCACTGCGGAATTGAAGTTGGAAGCCC
 TCACCCAGCGTCTGGAGCGAGAGATGGATGCTCACCCAAAGCGGACTACTTCGGTTCTGTGTGAAATG
 CAGCAAAGGGGTGTTGGAGCTGGCCAGGCCTGTCAGGCCATGGGAGATCTCTACCAGATGCGTGCTTC
 ACCTGTGCAGCCTGCAGCAGGAAGTTAAGAGGAAAGGCCTTCTATTTGTCAATGGCAAAGTATTCTGTG
 AGGAAGACTTTTTGTATTCTGGCTTTCAGCAGTCTGCGGACAGGTGTTTTCTTTGGACACCTGATCAT
 GGACATGATCCTACAGGCCCTAGGAAGTCTATCACCCGGCTGTTCCGCTGCGTCATCTGTAATGAA
 TGTCTGGATGGGGTTCCTTTCACCGTGGACTCTGAGAACAGATCTACTGTGTCGAGACTATCACAAGG
 TGCTGGCCCCCAAGTGTGCAGCCTGTGGCCTTCCATCCTTCCGCCAGAGGGTTCAGACGAGACCATCAG
 AGTTGTGTCCATGGACAGGGACTACCACGTAGAGTGCTACCAGTGTGAGGACTGTGGTCTGGAGCTCAAT
 GACGAGGACGGTCACCGCTGCTATCCACTGGAGGACCCTGTTCTGTCACTCCTGCCATGTCAAGAGAC
 TGGAGAAAGGACCCTCACCTGCACCCCTCCACCAGCACCATT**C**TA**G**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_013860

Insert Size: 2007 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_013860.2, NP_038888.2</u>
RefSeq Size:	4864 bp
RefSeq ORF:	2007 bp
Locus ID:	29806
UniProt ID:	<u>Q9QXD8</u>
Cytogenetics:	9 74.0 cM
Gene Summary:	Adapter or scaffold protein which participates in the assembly of numerous protein complexes and is involved in several cellular processes such as cell fate determination, cytoskeletal organization, repression of gene transcription, cell-cell adhesion, cell differentiation, proliferation and migration. Positively regulates microRNA (miRNA)-mediated gene silencing and is essential for P-body formation and integrity. Acts as a hypoxic regulator by bridging an association between the prolyl hydroxylases and VHL enabling efficient degradation of HIF1A. Acts as a transcriptional corepressor for SNAI1- and SNAI2/SLUG-dependent repression of E-cadherin transcription. Negatively regulates the Hippo signaling pathway and antagonizes phosphorylation of YAP1. Inhibits E2F-mediated transcription, and suppresses the expression of the majority of genes with E2F1-responsive elements. Regulates osteoblast development, function, differentiation and stress osteoclastogenesis. Enhances the ability of TRAF6 to activate adapter protein complex 1 (AP-1) and negatively regulates the canonical Wnt receptor signaling pathway in osteoblasts. May act as a tumor suppressor by inhibiting cell proliferation.[UniProtKB/Swiss-Prot Function]