

Product datasheet for **MC220489**

Samd4b (NM_175021) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	Samd4b
Synonyms:	BC042901
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



Fully Sequenced ORF:

>MC220489 representing NM_175021

Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGATGTTCCGAGACCAAGTGGGCATCCTCGCTAGCTGGTTCAAGGGCTGGAATGAGTGTGAACAGACAG
TGGCCCTGTGTCCCTTCTGAAGCGGGTCACCCGTACCCAAGCCCGCTTTCTGCAGCTCTGCCTGGAGCA
CTCGCTGGCGGACTGCAATGATATCCACCTGCTGGAGTCAGAGGCCAACAGTGTGCCATCGTCAGCCAG
TGGCAACAGGAGTCCAAGAGAAGGTGGTGTCCCTCTGTGTCCCATCTGCCCTGCTCCAGCCAGGCA
ACACGGAGGCCAAGTCCGAGTACATGAGACTGCTGCAGAAAGTGTGGCCTACTCCATCGAGAGCAATGC
CTTCATTGAGGAGAGCCGCAACTCCTCTCTATGCCCTCATCCACCCAGCCACCACACTGGAGGACCGC
AATGCGCTAGCCCTCTGGCTGAGCCACCTGGAAGAGCGGCTAGCTAGTGGTTTCCGCACCCGCCAGAGC
CCTCTACCATTCACGCCAGGGCTCAGATGAATGGGGTGGACCTGCAGAGCTGGCCCTGGAGAGGCAGG
GCCAGGCTGGCAGGACAAGCCACCCAGGGAATGGACACGTACCCCTCCACCCATCCAGCTCTGTGCCG
CCGGCCATCAACAGTATTGGGAGCAATGCAAAACACAGGTCTCCCTGCCAAATCCACCCTAGCCCACTGA
AGCGCTCCATGTCACTCATCCCTACGAGCCACAGGCCCTGGTGAAGTGGCCGAGTCCAGAGGAGCTTGG
GGCCCGGGCTGTCTTACCACGCCCCGACCACGCCCTCTCGCCACAGAGCAGTGTGGCTCTCTGTGGC
AGTGAGCAGACGGAGGAGCAGGGCTCCAGCCGGAACACTTTCCAGGAGGATGGCAGCGGGATGAAAGATG
TGCCCTCATGGCTCAAGAGCCTCCGCTTGCAAGATGCTGCTCTGTTCTCACAGATGAGCTACGAAGA
GATGATGACCTCACTGAGCAGCACCTAGAGTCGAGAATGTACCAAAGGTGCCCGGCACAGATAGCC
CTGAGCATCCAAAGCTGCGGGAGAGACAGAGCGTCTCAAGTCCCTAGAGAAGGATGTGCTGGAAGGTG
GCAATCTCGGAATGCTCTGCAGGAAGTGCAGCAGATCATCACTCTATCAAGGCCTACAGTGTCTCT
TCAGGCCACCCCTACTGCCAAGGATGAGGGTCGGGGGAGCCGCTGCTGCCAGGTGCTGAGCCTCCCTC
ACCCACCTGGATCAGACAAGGGCACTGAGGCCAAGGACCTCCAGCTGCAGAGAACTATCTCTCCAC
CAGCTCCAGCCCCCTCTGATAGCAGTGAGCCAGCCAGCTCCCGTCGCCGACGGAGACATCCCAGCCA
GTTTACACGGGTGATGGCAAGTATGCACCCAGCTGCTGGTGTCCCGACCAGACGAGGAGAACATCACC
AGTTACCTCCAGCTCATCGAAAAGTGCCTGACTCATGAGGCTTTCACAGAGACACAGAAGAAGCGTTGC
TGTCCTGGAAGCAGCAAGTGTGAACTCCTCCGACATTCCTCGAAAAGCTGCACTGGACATGCAGAG
CTACCGGCAGCAGAAGGGCTGGGCATTTGGCTCCAACCTCGTCCCATAGCTGGCTCTGTGGGGATGGGT
GTTGCCCGGGACCCAGCGCAGTTCCTCAATGCCTCTCGGGCATTACCACCTGGCAGGATGGGCCTCC
TGAGCCCTTCAGGCATTGGGGTGTCTCTCTCGACATGCCCTTACCAGCCCCAGTCTTGGGGGTGAGGG
CCGACAGAACCTGTGTTTGCATCCTGGAGGCAGCAACAGCATGCCAGCCAGAGCCGACGCTCTGTG
CAGCGCACTCACTCGTCCCAGTTCCTCGTCAACCCAGGCTATTCTCATGTTCCCTCCAGACTGCCCGG
TACCTGGGCCTGACCTGGAGATCAATCCTACCTGGAGTCTCTGTGTCTGAGCATGACAGAACACGCCTT
GGGTGATGGGACAGACAAAACCTCCACCATCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:

SgfI-MluI

ACCN:

NM_175021

Insert Size:

2064 bp

OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery. The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_175021.3 , NP_778186.2
RefSeq Size:	4391 bp
RefSeq ORF:	2064 bp
Locus ID:	233033
Cytogenetics:	7 B1