

Product datasheet for **MC221435**

Ubtf (NM_011551) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ubtf (NM_011551) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Ubtf
Synonyms:	A930005G04Rik; NOR-90; Tcfubf; UBF; UBF-1; UBF1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >MC221435 representing NM_011551
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGAACGGAGAAGCGGACTGTCCCACAGACCTGGAATGGCCGCCCCAAAGGCCAAGACCCTGGTCCC
 AGGAAGATATGCTGACTCTGCTGGAATGCATGAAGAACAACCTCCCATCCAACGACAGCTCCAAGTTCAA
 AACACAGAGTCCACATGGATTGGGAAAAAGTTGCGTTTAAAGATTTTTAGGAGATATGTGCAAGCTC
 AAATGGGTGGAGATTTCTAACGAGGTAAGGAAGTCCGTACATTGACAGAATTGATCCTTGATGCTCAGG
 AACATGTTAAAAACCTTACAAAGGCCAAAAAGCTCAAGAAACACCCGGATTTTCCAAAGAAACCTCTCAC
 CCCTTACTTCGCTTCTTATGGAGAAGCGGGCCAAAGTACGCGAAGCTCCACCCGGAGATGAGCAACCTG
 GACCTGACTAAGATCCTGTCTAAGAAATACAAGGAAGTCCAGAGAAGAAGAAGATGAAATATATTCAGG
 ACTTCCAGAGGGAGAAACAGGAGTTCGAGCGAAACCTGGCCCGATTCCAGGGAGGATCACCTGACCTTAT
 CCAGAATGCCAAGAAGTCGGACATCCCCGAGAAACCCAAACTCCCAGCAACTGTGGTACACCCACGAG
 AAGAAGGTGTATCTCAAAGTGGCGCCGGATGCCACTACGAAGGAGGTGAAGGACTCCCTGGGGAAGCAGT
 GGTCTCAGCTCTCGGACAAAAAGAGGCTGAAATGGATTATAAGGCCCTGGAGCAGCGGAAGGAGTACGA
 GGAGATTATGCGGGACTATATCCAGAAGCACCTGAGCTGAACATCAGTGAGGAAGGCATCACCAAGTCC
 ACCCTCACCAAGGCCGAGCGCCAGCTCAAGGACAAGTTCGATGGGCGACCTACCAAGCCGCTCCGAACA
 GCTACTCCCTGTACTGTGCGGAGCTCATGGCCAACATGAAGGATGTGCCAGCACAGAGCGCATGGTGCT
 ATGCAGCCAGCAATGGAAGCTGCTCTCTCAGAAGGAGAAGGATGCCTATCACAAGAAATGTGACCAGAAA
 AAGAAAGATTATGAGGTAGAGCTGCTGCGGTTCCCTTGAGAGCTTGCCGGAGGAAGAGCAGCAGCGGGTTC
 TGGGGGAGGAGAAGATGCTGAACATCAATAAGAAGCAAACACCAGTCCGGCCCTCAAGAAGCCTTCACA
 GGAAGGTGGCAAGGGTGGCTCCGAGAAGCCAAAGCGGCCTGTGTCTGCTATGTTTCACTTCTCTGAGGAG
 AAGCGAAGGCAGCTACAGGAGGAACGACCTGAGCTCTCCGAAAGCGAACTGACCCGCCTTCTGGCCCGCA
 TGTGGAACGACTTGTCCGAGAAGAAGAAGGCTAAATACAAGGCCCGGGAGGCTGCGCTGAAGGCGCAGTC
 TGAGAGGAAGCCGGGCGGGGAGCGTGAAGATAGGGGCAAGCTGCCGGAGTCGCCAAGAGGGCTGAGGAG
 ATCTGGCAGCAGAGTGTATCGGAGACTATCTGGCCCGCTTCAAGAATGACCGGGTGAAGCCTTGAAG
 CCATGGAGATGACTTGAACAACATGGAGAAGAAGGAGAAGCTGATGTGGATTAAGAAGGCAGCAGAAGA
 CCAAAAAACGATATGAGAGAGAGTTAAGTGAGATGCGGGCCCTCCAGCTGCTACAAACTCTTCAAGAAG
 ATGAAGTTCAGGGAGAACCAAGAAACCTCCTATGAACGGTTACCAGAAGTTCTCCAGGAGCTGCTGT
 CCAATGGGGAGCTAAATCACCTGCGCTGAAGGAGCGCATGGTGGAGATTGGCAGCCGCTGGCAGCGCAT
 CTCCCAGAGCCAGAAGGAGCACTATAAGAAGCTGGCGGAGGAGCAACAGAGGCAGTACAAGGTGCACTTG
 GACCTCTGGGTCAAGAGCCTGTCTCCCCAGGACCGCGCAGCATACAAAGAATACATCTCCAATAAACGTA
 AGAACATGACCAAGCTCCGAGGCCAAACCCCAAGTCTAGCCGGACCACCTGCAGTCCAAGTCCGGAGTC
 CGAGGAGGATGACGATGAGGAGGAGGAGGACGACGAGGAGGAGGAGGAAGAGGAAAGATGATGAGAACGGT
 GACTCTTCTGAGGATGGCGGGGACTCTTCTGAGTCCAGCAGTGAAGATGAAAGCGAGGATGGGGATGAGA
 ATGATGATGACGACGACGATGAAGATGACGAAGATGACGACGATGAGGATGAAGACAACGAGTCTGAGGG
 CAGTAGCTCCAGCTCTTCACTCAGGGGACTTTCGGATTCTGACTCCAAC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja2604_c02.zip

Restriction Sites: SgfI-MluI

ACCN: NM_011551

Insert Size: 2295 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_011551.6](#), [NP_035681.2](#)

RefSeq Size: 4704 bp

RefSeq ORF: 2295 bp

Locus ID: 21429

Cytogenetics: 11 D