

Product datasheet for **MC222029**

Atg9a (NM_001003917) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Atg9a (NM_001003917) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Atg9a
Synonyms:	Apg9l1; Atg9; Atg9l1; AU019532
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >MC222029 representing NM_001003917
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCACAGTTTGACACTGAATACCAGCGCCTAGAGGCCTCTATAGCGATTGCCCCCGGGGAAGAGG
 ACCTGTTGGTGCATGTGGCCGAGGGGAGCAAATCACCTTGGCACCACATTGAAAACCTTGACCTTCTT
 CTCTCGAGTTTATAATCTCCACCAGAAGAATGGCTTACATGTATGCTCATTGGAGAGATGTTTGAACCT
 ATGCAGTTCCTCTTTGTGGTTGCCTTACCACCTTCTGGTTAGCTGTGTGGACTATGACATCCTCTTTG
 CCAACAAGATGGTGAACCACAGTCTCCATCCTACTGAGCCTGTCAAGGTCCTGCCAGATGCCTTTTT
 GCCTGCCCAAGTCTGTAGTGCCAGGATTCAGGAAAATGGCTCTTATCACCATCCTGGTCATTGCTGGT
 GTCTTCTGGATCCACCGCTTATCAAGTTTATCTATAACATTTGCTGCTATTGGGAGATCCACTCCTTCT
 ACCTACATGCTCTTCGGATCCCAATGTCTGCCCTTCCATACTGCACATGGCAGGAAGTTCAGGCCGGAT
 TGTGCAGACGCAGAAAGAGCATCAGATCTGCATCCACAAGCGTGAGCTGACAGAGTTGGACATCTACCAT
 CGCATCCTACGTTTCCAGAACTACATGGTGGCCCTGGTGAACAAATCCCTCCTGCCTCTGCGCTTCCGTC
 TGCCCCGGCCTCGGAGAGGTTGTCTTCTTCAACCGTGGCCTCAAGTACAACCTTTGAGCTCATCCTCTTCTG
 GGGACCCGGCTCTCTGTTTCTCAATGAGTGGAGCCTCAAGGCCGAGTACAAACGTGGAGGGCAACGGCTA
 GAGCTGGCCCAGCGTCTCAGCAACCGCATCTTGTGGATTGGCATCGCCAACTTCTGCTGTGTCCCTCA
 TCCTCATCTGGCAGATCCTCTATGCCTTCTTACGCTATGCCGAGGTGCTGAAGAGAGAGCCGGGGCCCT
 GGGAGCGCGTTGCTGGTCACTCTATGGCCGTTGTTACCTCCGCCACTTCAATGAGCTGGAGCATGAGCTG
 CAGTCCCGCCTCAACCGAGGCTACAAGCCCGCCTCCAAGTACATGAATTGCTTCTTGTACCAGGCTGCTGA
 CTCTGTGGCCAAGAATGGTGCCTTCTTCGCTGGCTCTATCCTGGCTGTGCTTATTGCCCTCACCATCTA
 TGATGAAGATGTGTTAGCTGTGGAACACAGTCTCACCACGGTCAACCTCCTGGGAGTACGGTGACTGTG
 TGCAGGTCCTTATCCCAGACCAGCACATGGTGTCTGCCCGAGCAGCTGCTCCGAGTATTCTTGCAC
 ACATCCACTACATGCCTGACCACTGGCAGGGTAAATGCCACCGCTCGCAGACCCGGGACGAGTTTGCCCA
 GCTCTTCCAGTACAAGGCAGTGTTCATCTTGGAGGAGTTGCTGAGTCCCATCGTCACACCCCTCATTCTC
 ATCTTCTGCCTCCGCCCTCGGGCCTGGAGATCATAGACTTCTTCCGCAACTTACGGTCGAGGTCGTGG
 GTGTGGGAGACACCTGCTCCTTGTCTCAGATGGACGTTCCGACAGATGGCCATCCTCAGTGGCTGTCTGG
 AGGGCAGACAGAGGCTCAGTGTACCAGCAAGCCGAGGACGGGAAGACCGAGTTGTCGCTCATGCACTTT
 GCCATACCAATCCCGGCTGGCAGCCCCCTCGTGAGAGCACAGCTTTCCTGGGCTTCTCAAGGAGCAGG
 TGCAGCGAGATGGAGCAGCTGCTGGCCTGGCCAGGGTGGCCTGCTCCCCGAGAATGCCCTTTCACATC
 CATCCAGTCTTACAGTCTGAGTCCGAGCCACTGAGCCTTATTGCAAATGTGGTAGCAGGCTCATCCTGC
 CGAGGACCCTCACTGTCCAGAGACCTGCAGGGCTCCAGGCACAGGGCTGATGTTGCTTCTGCCCTTCGAT
 CCTTCTCCCCTCTGCAGCCTGGAGCGGCCCTCAAGGCCGGGTCCCAGTACCATGACAGGCTCTGGAGT
 GGATGCCAGGACAGCCAGCTCTGGGAGTAGCGTGTGGGAAGGACAGCTGCAGAGCCTGGTGTGTCCGAA
 TACGCGTCCACCGAGATGAGCCTGCACGCCCTGTACATGCACCAGCTCCACAAGCAGCAGACCCAGGCTG
 AGCCCGAGCGGCATGTGTGGCACCGCCGGGAGAGTGTGAGAGTGGAGAGAGTGGCCCTGAAGAGGGGGG
 AGAGGGTGGCCGGCCCCCAACCCATCCCCCGCTCGGCCAGCTATCCCTGTGCTACACCCCGCCTGGA
 GCACCTGAGACCACCGCCTGCATGGGGCTTCCAGAGGCGCTACGGGGGCATCACAGATCCTGGCACAG
 TGCCCCGTGGCCCTCTCACTTCTCAGGCTGCCCTGGGAGGATGGGAGAAAGATGGCCAGCCAGCATC
 AAGGCACCCAGAGCCGGTCCAGAGGAGGGCTCAGAGGATGAACTCCCCCTCAAGTGCACAAGGT**ATAG**

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-RsrII
ACCN: NM_001003917
Insert Size: 2520 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001003917.4</u> , <u>NP_001003917.2</u>
RefSeq Size:	3856 bp
RefSeq ORF:	2520 bp
Locus ID:	245860
UniProt ID:	<u>Q68FE2</u>
Cytogenetics:	1 C4
Gene Summary:	<p>Involved in autophagy and cytoplasm to vacuole transport (Cvt) vesicle formation. Plays a key role in the organization of the preautophagosomal structure/phagophore assembly site (PAS), the nucleating site for formation of the sequestering vesicle. Cycles between a juxta-nuclear trans-Golgi network compartment and late endosomes. Nutrient starvation induces accumulation on autophagosomes. Starvation-dependent trafficking requires ULK1, ATG13 and SUPT20H (By similarity). Required for carbonyl cyanide m-chlorophenylhydrazone (CCCP)-induced ATG8 family proteins lipidation, a key autophagy step.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1 and 2 both encode the same isoform (a).</p>