

Product datasheet for **MC223121**

Zcchc6 (BC023880) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zcchc6 (BC023880) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Zcchc6
Synonyms:	MGC38585, 6030448M23Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC023880 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGTTTTCTTCCTTCAGCAAAGAAAAGAACCTCTCTTGCCAGTTTATCTAGGATCATGGATTGAAG
AATTCTCATTAAATAAATTAGGGAATTTACGCCTTAAAGATGTTGAGAAAGACTCTGTGGTCTGGGAATA
CACTGACAACAGCACAGGGGACACCAGCTCAGCAAAAGAGGAGGCACCCAAAGAGACTGCAGCAAAGAAG
GGACAGGTGCCATTGACATTTAATAAAGCATCAGCCTTCAGTGCCAGTCGGGCAACTCTGGGTGGAAT
TGCTGCGGTTCTATGCTTTAGAGTTTAATCTGGCTGATTTAGTGATCAGTATTCGTGTCAAAGAATTAAT
ATCACGGGAATCAAAGGACTGGCCCAAAAAGCGTATCGCAATTGAAGACCCCTATTCTGTAAAAAGAAAT
GTGGCAAGGACCCCTGAATAACCAAGCCTGTGTTGAGTATATACTCCATTGTTTGAGGACAACATATAAAT
ATTTGCGCTCTCCACACAAAGTTACAAAACCAACCTCACAAAGCCTCCTAGCCAGTCACCTGTGTTTC
AGACCCCTACAGAGAAGCAAAAATGGTGGCCCAAGCAACAAATATTGATAAACTGGGAAAT
GCAGCTGTGGCTCAAGATCCTGGTGTCAAACCTCAGGTGACTGCAGAGCACAGCTAGTTACTCTTAAGA
ATACCACTGAGGAGGTGGGAAGCCAGCAAGGAGAAGACAGGAGGTGTGCATATCCAGCCACCAAGGA
GAGCTCTGGCTGTGTCCAGGCAGAGGTGAGTTGTGAGGGTTTGGAGGATGCCACAGCAGAGCTTCGGAA
ACAGGAAGTGACAACGAGGAAGTCAGAAGGAAGACCAAGCATCCTTTGAGCACTGATGACCAAGGTTTAA
GCAGCAGCAAGCACCCAGAACTGCAGAACTGTGGCAGCCTGTGTGGTCTGCAGGCAGACAACACTCTGGA
GTTAGTAGCTGAAGAGTGAATAGCTGTGCCTCTTTGGACAACAAGGCTGAGGTCAATGAAGAGAGGATA
GAAGGCGCTGAGGAAGTGAAGAAGCAGCAGCTCTGAGTTGCTTTTCCCCTTCAGTACAGAGTCGCACAT
CTGCAGCCATGCACTTTGATGATGAGGAGGAGGAGGAGGAAGAAGAGGAGGAAGAGGAGCCAGGCTCAG
CATTAACCTTACAGAAGATGAGGAGGGGTTGCTAATGAGCATCAGGTGGACAGCAGGTATGCTGGGTCA
GGTGAGGAGGATGCCCTCTCTGAAGAAGATGACTTGGCAGAGCCTGCTAAGGGTGAAGATACAGGAGAAT
GTGGAGAAAACGTGGGTGGGACTCTGCTCATAGACCTAAATAGAATAACACTTAAAGAAGAAAGCTTTCC
TGAGGAAGATTTACCTGGGACCAGTCTGAGTTCTTCTACGAATTTAGAAAAGTACCTTCACCAAAGGC



[View online »](#)

AAGTCCCCACAGTAGTGTGCAGCTTGTGTAAGCGAGAAGGTCACCTAAAGAAGGACTGTCCTGAGGACT
TCAAGAGGGTCCAGCTGGAGCCGTTGCCACCCCTGACACCCAAGTTCTCAAACATCTTAGATCAAGTCTG
CGTCCAGTGTTACAAGGATTTTCTCCTCACTATTGTAGAAGATCAGGCTCGTGAACATATCCGGCAAAAT
CTAGAAAGTTTCATAAAACAGGATTTTCCAGGAACTAAATTGAGCTTGTTGGCTCCTCCAAAAATGGAT
TTGGGTTCAAACAGAGTGACCTTGACGTCTGTATGACGATTACCGACATGAGACTGCTGAGGGGCTGGA
CTGTGTAAGAACTATTGAGGAGCTGGCTCGAGTCTCAGAAAGCACTCAGGGTTGAGGAATATCCTTCCT
ATCACAACAGCCAAGGTGCCATCGTGAAGTTCTCCATCTAAGAAGTGGTCTGGAAGTGGACATCAGCT
TGTATAACACACTGGCGCTTCATAACACAAGGCTCTGTCTGCATACTCAGCCATTGACCCAGAGTGAA
GTATCTGTGCTACACCATGAAAGTGTTCAAAAGATGTGTGACATCGGTGACGCTCCAGAGGCAGCTTG
TCATCCTACGCTTACACTCTTATGGTGCTATACTTCTCCAGCAAAGGTCTCCACCTGTCATCCCGTGC
TTCAAGAGATATACAAAGGTGAAAAGAAACCAGAAATACTTGTGATGGCTGGAATATTTACTTTTCGA
TCAAATAAATGAACTGCCACTTGTGGCCAGAATATGGGAAAAATACAGAACCTGTGGGCAGCTGTGG
CTGGGACTTCTGCGCTTCTACACTGAGGAGTTTGACTTTAAGGAGCAGTGATCAGCATCAGAAGAAAGA
GTCTGCTTACAACCTTTAAGAAACAGTGGACCTCAAAATACATTGTTATTGAAGATCCCTTCGATTTGAA
TCACAATCTTGGTGCTGGATTATCGAGGAAAATGACAAATTTTATAATGAAAGCTTTTATCAATGGTAGA
AGAGTATTTGGAATTCAGTCAAAGGATTTCAAAGGACAACCCCTCAAACTGGCATACTTTTTTGATC
CGGATGTCTTAAGTGAAGGAGAACTGGCCCCAAATGACAGGTGCTGCCGAATTTGTGGGAAAATCGGGCA
CTTCATGAAGGACTGTCCATGAGGAGGAAAGTGAAGGCGACGGAGAGACCAGGAAGATACCCCAAACCA
AGGTACTCTGAGAGCAAAGAAAAAGAAAGCAAAGAGGACAAAGAGATTGAGAACAAGTACACAGAAAAGG
AAGTGTCAACAAAAGAGATAAGCTCACACCGTGTGCAGCTGCGAAAGCCAAGCCAGTGAGGGCAGCTGT
CGACCTGGGCGGGGAGAGCTCCTCAGGACGCCACAGAAAAGTGAAGAGACAGGACGACAGAGACTCG
AGAGAAAAGCGTTGTTTCATCTGTGGAAGAGAAGGACACATTAAGGAATGCCACAGTTCAAAGGCT
CTCCAGGTAGCCTTTCCAGTAAATATGACTCAGGGGAGAGCCTCAGTGAAGAGGACCCAGCAGGAATC
ATGA

AGCGGACCGACGCTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATA
TCCTGGATTACAAGGATGACGACGATAAGGTTAA

Restriction Sites:	SgfI-RsrII
ACCN:	BC023880
Insert Size:	3084 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC023880](#), [AAH23880](#)

RefSeq Size: 4215 bp

RefSeq ORF: 3083 bp

Locus ID: 214290

Cytogenetics: 13 B2

Gene Summary: Uridyltransferase that mediates the terminal uridylation of mRNAs with short (less than 25 nucleotides) poly(A) tails, hence facilitating global mRNA decay (PubMed:28792939). Essential for both oocyte maturation and fertility. Through 3' terminal uridylation of mRNA, sculpts, with TUT7, the maternal transcriptome by eliminating transcripts during oocyte growth (PubMed:28792939). Involved in microRNA (miRNA)-induced gene silencing through uridylation of deadenylated miRNA targets. Also acts as a suppressor of miRNA biogenesis by mediating the terminal uridylation of miRNA precursors, including that of let-7 (pre-let-7) (PubMed:22898984). Uridylated pre-let-7 RNA is not processed by Dicer and undergo degradation. Pre-let-7 uridylation is strongly enhanced in the presence of LIN28A. Due to functional redundancy between ZCCHC6 and ZCCHC11, the identification of the specific role of each of these proteins is difficult (By similarity) (PubMed:22898984). Involved in microRNA (miRNA)-induced gene silencing through uridylation of deadenylated miRNA targets (By similarity). Also functions as an integral regulator of microRNA biogenesis using 3 different uridylation mechanisms (By similarity). Acts as a suppressor of miRNA biogenesis by mediating the terminal uridylation of some miRNA precursors, including that of let-7 (pre-let-7). Uridylated pre-let-7 RNA is not processed by Dicer and undergo degradation. Pre-let-7 oligouridylation is strongly enhanced in the presence of LIN28A (PubMed:22898984). In the absence of LIN28A, TUT7 and TUT4 monouridylate group II pre-miRNAs, which includes most of pre-let7 members, that shapes an optimal 3' end overhang for efficient processing (By similarity). Add oligo-U tails to truncated pre-miRNAs with a 5' overhang which may promote rapid degradation of non-functional pre-miRNA species (By similarity). Does not play a role in replication-dependent histone mRNA degradation (By similarity). Due to functional redundancy between TUT4 and TUT7, the identification of the specific role of each of these proteins is difficult (PubMed:28792939, PubMed:22898984). TUT4 and TUT7 restrict retrotransposition of long interspersed element-1 (LINE-1) in cooperation with MOV10 counteracting the RNA chaperone activity of L1RE1. TUT7 uridylates LINE-1 mRNAs in the cytoplasm which inhibits initiation of reverse transcription once in the nucleus, whereas uridylation by TUT4 destabilizes mRNAs in cytoplasmic ribonucleoprotein granules (By similarity).[UniProtKB/Swiss-Prot Function]