

Product datasheet for MC223439

Xpo1 (NM_001035226) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Xpo1 (NM_001035226) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Xpo1
Synonyms:	AA420417; Crm1; Exp1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>MC223439 representing NM_001035226 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGCCAGCAATTATGACAATGTTAGCAGACCATGCAGCTCGTCAGCTGCTTGATTTAGCCAAAACTGG
ATATCAACCTTTTAGACAACGTGGTGAATTGCTTATACCATGGAGAAGGAGCCAGCAAAGAATGGCTCA
AGAAGTACTGACACATTTAAAGGAGCATCTGATGCTTGGACACGGGTCGACACAATATTGGAGTTTTCG
CAGAACATGAACACGAAATACTATGGACTACAGATTTTGGAAAAATGTGATCAAAAACAGGTGGAAGATCC
TTCCAAGGAACAGTGTGAAGGAATAAAAAATATGTTGTTGGTCTCATTATCAAGACATCCTCTGACCC
AACTTGTGTAGAGAAAAGAAAAGGTGATATTGGGAAATGAATATGATTCTTGTTCAGATACTTAAGCAA
GAGTGGCCTAAGCATTGGCCAACTTTTATCAGTGACATTGTTGGAGCGAGTAGGACCAGTGAAAGTCTAT
GTCAAAAATAATGTTGATTCTTAAACTCTTGAGTGAAGAAGTGTGATTTCTCTAGTGGACAGATAAC
CCAAGTTAAAGCAAAGCATTAAAAGACAGCATGTGCAATGAGTTTTACAAAATATTTAGCTTTGTGTCAG
TTTGTGATGGAAAATCCCAAAATGCTCCCTTGTTCATGCTACATTAGAAACACTGCTGAGATTTCTGA
ACTGGATTCCACTGGGATACATCTTTGAGACCAAGTTAATCAGCACATTAATTTACAAGTTCCTGAAATG
TCCAATGTTCCGAAATGTGCTCTAAAGTGCCTCACTGAGATTGCTGGTGTGAGTGAAGCCAATATGAG
GAACAATTTGAAACTCTGTTTACACTTACAATGATGCAGCTAAAACAGATGCTTCTTTAAACACTAATA
TCCGACTTGCATACTCGAATGGGAAAGATGATGAACAGAACTTTATTCAAAATCTCAGCCTGTTTCTCTG
TACCTTTCTAAAAGAGCACGGTCAACTTCTAGAAAAACGACTGAATCTCAGGGAAGCACTCATGGAGGCT
CTTCATTATATGTTGTTGGTGTCTGAAGTGAAGAAACTGAAATTTTTAAAATCTGCCTTGAGTACTGGA
ATCATTGGCAGCTGAGCTCTACAGAGAGAGTCCATTTTCTACATCTGCTTCTCCCTTGTATCTGGAAG
TCAACACTTTGATTTCTCCTAGGAGACAGCTGTAAGTACCCTGTTATCAAAGTCCGTTTATTGATG
GTTAGTGCATGGCTAAACCAGAGGAGTACTGGTTGTAGAAAAATGATCAGGGAGAAGTTGTAAGAGAGT
TCATGAAGGATACTGATTCCATAAATTTATATAAGAATATGAGAGAAACATTAGTTTATCTTACTCATCT
GGATTATGTAGATACGGAATAATAATGACTAAGAACTTCAAAATCAAGTGAATGGTACAGAATGGTCA
TGAAGAATTTGAATACGTTGTGTTGGCAATAGGCTCCATCAGTGGCGCAATGCATGAAGAGGATGAGA



[View online »](#)

```

AGCGGTTTCTTGTTACTGTTATAAAGGATCTTTTAGGATTATGTGAGCAGAAAAGAGGCAAAGATAATAA
AGCTATTATTGCATCGAATATCATGTACATAGTAGGTCAATACCCAAGATTTTTGAGAGCACACTGGAAA
TTTTTGAAGACAGTAGTAAACAAGCTATTCGAATTCATGCATGAGACCCATGATGGAGTCCAGGATATGG
CTTGTGACACATTTATTTAAAATAGCTCAAAAATGCCGTAGGCATTTTGTTCAGGTTCAAGTTGGAGAAGT
AATGCCGTTTCATTGATGAGATTTTAAACAACATAAATACTATAATTTGTGATCTTCAGCCTCAGCAGGTT
CATACATTTTATGAAGCTGTGGGTACATGATTTGGTGCACAGACAGACCAGACAGTTCAAGAGCATTGGA
TAGAAAAATACATGCTACTCCCTAATCAAGTGTGGGACAGTATAATTCAGCAGGCTACCAAAAATGTGGA
TATCCTGAAAGATCCTGAAACAGTCAAGCAGCTGGGTAGCATATTTAAAAACAAATGTGAGAGCCTGCAAA
GCTGTTGGTCATCCCTTTGTAATTCAGCTTGGCAGAATTTATTTAGACATGCTTAATGTATATAAGTGCC
TCAGTGAAAAATTTCTGCAGCTATTCAAGCTAATGGTGAGATGGTTACAAAGCAACCATTGATTAGAAG
TATGCGGACAGTAAAAAGGGAACTTTAAAGTTAATATCTGGTTGGTGAGCCGTTCCAATGATCCTCAG
ATGGTAGTGAAAACTTTGTTCCCTTTGTTGGATGCAGTTCTTATTGATTACCAGAGAAATGTACCAG
CTGCTAGAGAACCAGAAGTGCTTAGTACTATGGCCATTATTGTCAACAAATTAGGTGGACATATAACAGC
TGAATAACCTCAAATATTCGATGCAGTTTTTGTGAGTGCACATTGAATATGATAAATAAGGATTTTGAAGG
TATCCAGAACACAGAACAACTTTTTCTTATTACTTCAAGCTGTCAATTCTCACTGTTTCCAGCATTCC
TAGCTATTTCCACCTGCACAGTTCAGCTGGTTTTGGATTCTATTTTGGGCTTTTAAACATACTATGAG
GAATGTTGCAGATACAGGCTTACAGATACTTTTTACACTTTTACAAAATGTTGCACAAGAGGAAGCTGCC
GCACAGAGTTTTTATCAAATTTTGTGATATTCTTCAGCATATTTCTCTGTTGTAACAGACACTT
CACATACTGCAGGTCTAACAATGCATGCATCAATTCTTGCAATATGTTTAAATTTGGTTGAAGAAGGAAA
AATAAGTACACCACTAAATCCTGGAAATCCAGTTAACAACCAATGTTTATTCAAGACTATGTGGCTAAT
CTCCTTAAATCTGCATCCCTCATCTACAAGATGCTCAAGTAAAGCTCTTGTAAACAGGACTTTTCAGCT
TAAATCAGGATATTCTGCTTTCAAAGAACACCTAAGGGATTTCTTAGTTCAAATAAAGGAGTTTGCAGG
TGAAGATACATCGGATCTATTTTGAAGAAAGAGAAACAGCCCTTCGACAGGCTCAGGAGGAGAAACAT
AACTTCAAATGTCTGTACCTGGGATTCTTAATCCACACGAAATTCAGAAAGAAATGTGTGATTA

```

```

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

```

Restriction Sites:

SgfI-MluI

ACCN:

NM_001035226

Insert Size:

3216 bp

OTI Disclaimer:

Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001035226.1](#), [NP_001030303.1](#)

RefSeq Size: 4906 bp

RefSeq ORF: 3216 bp

Locus ID: 103573

UniProt ID: [Q6P5F9](#)

Cytogenetics: 11 A3.2

Gene Summary: Mediates the nuclear export of cellular proteins (cargos) bearing a leucine-rich nuclear export signal (NES) and of RNAs. In the nucleus, in association with RANBP3, binds cooperatively to the NES on its target protein and to the GTPase Ran in its active GTP-bound form. Docking of this complex to the nuclear pore complex (NPC) is mediated through binding to nucleoporins. Upon transit of a nuclear export complex into the cytoplasm, disassembling of the complex and hydrolysis of Ran-GTP to Ran-GDP (induced by RANBP1 and RANGAP1, respectively) cause release of the cargo from the export receptor. The directionality of nuclear export is thought to be conferred by an asymmetric distribution of the GTP- and GDP-bound forms of Ran between the cytoplasm and nucleus. Involved in U3 snoRNA transport from Cajal bodies to nucleoli. Binds to late precursor U3 snoRNA bearing a TMG cap (By similarity).
[UniProtKB/Swiss-Prot Function]