

Product datasheet for MC223644

Tek (NM_013690) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tek (NM_013690) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Tek
Synonyms:	AA517024; Cd202b; Hyk; STK1; Tie-2; Tie2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC223644 representing NM_013690 Red=Cloning site Blue=ORF

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACTCTTTAGCCGGCTTAGTCTCTGTGGAGTCAGCTTGCTCCTTTATGGAGTAGTAGAAGCGCCA
TGGACCTGATCTTGATCAATTCCTACCTCTTGCTGTGATGCCGAAACATCCCTCACCTGCATTGCCTC
TGGGTGGCACCCCATGAGCCCATCACCATAGGAAGGGACTTTGAAGCCTAATGAACCAGCACCAAGAT
CCACTGGAGGTTACTCAAGATGTGACCAGAGAATGGGCGAAAAAGTTGTTTGAAGAGAGAAAAAGCCA
GTAAGATTAATGGTGCTATTTCTGTGAAGGTCGAGTTCGAGGACAGGCTATAAGGATACGGACCATGAA
GATGCGTCAACAAGCATCCTTCTACCTGCTACTTTAACTATGACCGTGGACAGGGGAGATAATGTGAAC
ATATCTTTCAAAAAGGTGTTAATTAAGAAGAAGATGCAGTGATTTACAAAAATGGCTCCTTCATCCACT
CAGTGCCCGGCATGAAGTACCTGATATTTTAGAAGTTCACTTGCCGCATGCTCAGCCCCAGGATGCTGG
TGTGTAAGTCTCAGAAAGTGGGGCCCGACTGTAGCCGTCTTGTACTACTTGCAAGAACAATGGAGTCTGCC
ATGAAGATACCGGGGAATGCATTTGCCCTCCTGGGTTTATGGGGAGAACATGTGAGAAAGCTTGTGAGCC
GCACACATTTGGCAGGACCTGTAAGAAAGGTGTAGTGGACCAGAAGGATGCAAGTCTTATGTGTCTGT
CTCCCAGACCTTACGGGTGTTCTGTGCCACAGGCTGGAGGGGTTGCAGTGAATGAAGCATGCCCAT
CTGGTTACTACGGACAGACTGTAAGCTCAGGTGCCACTGTACCAATGAAGAGATATGTGATCGGTTCCA
AGGATGCCTCTGCTCTCAAGGATGGCAAGGCTGCAAGTGTGAGAAAGAAGGCAGGCCAAGGATGACTCCA
CAGATAGAGGATTTGCCAGATCACATTGAAGTAAACAGTGGAAAAATTAACCCCATCTGCAAAGCCTCTG
GGTGGCCACTACCTACTAGTGAAGAAATGACCCTAGTGAAGCCAGATGGGACAGTCTCAACCAAAATGA
CTTCAACTATACAGATCGTTTCTCAGTGGCCATATTCAGTGTCAACCGAGTCTTACCTCTGACTCAGGA
GTCTGGGCTGCAAGTGTGAACACAGTGGCTGGGATGGTGGAAAAGCCTTTCAACATTTCCGTCAAAGTTC
TTCCAGAGCCCTGCACGCCCAAATGTGATTGACACTGGACATAACTTTGCTATCATCAATATCAGCTC
TGAGCCTTACTTTGGGGATGGACCCATCAAATCCAAGAAGCTTTTCTATAAACCTGTCAATCAGGCTGG



```

AAATACATTGAAGTGACGAATGAGATTTTCACTCTCAACTACTTGGAGCCGCGGACTGACTACGAGCTGT
GTGTGCAGCTGGCCCGTCTGGAGAGGGTGGAGAAGGGCATCCTGGGCTGTGAGACGATTTACAACAGC
GTCTATCGGACTCCCTCCTCCAAGAGGTCTCAGTCTCCTGCCAAAAAGCCAGACAGCTCTAAATTTGACT
TGGCAACCGATATTTACAAACTCAGAAGATGAATTTTATGTGGAAGTTCGAGAGGGCGATCCCTGCAAAACA
CAAGTGATCAGCAGAACATCAAAGTGCCTGGGAACCTGACCTCGGTGCTACTGAGCAACTAGTCCCCAG
GGAGCAGTACACAGTCCGAGCTAGAGTCAACACCAAGGCGCAGGGGGAGTGGAGTGAAGAAGTCAAGGGCC
TGGACCCCTTAGTGACATTCTCCCTCAACCAGAAAACATCAAGATCTCCAACATCACTGACTCCACAG
CTATGGTTTCTTGACAATAGTGGATGGCTATTGATTTCTCCATCATCATCCGGTATAAGGTTTCAGGG
CAAAAAATGAAGACCAGCACATTGATGTGAAGATCAAGAATGCTACCGTTACTCAGTACCAGCTCAAGGGC
CTAGAGCCAGAGACTACATACCATGTGGATATTTTGTGAGAACAACATAGGATCAAGCAACCCAGCCT
TTTCTCATGAACTGAGGACGCTTCCACATTCCCCAGCCTCTGCAGACCTCGGAGGGGGAAAGATGCTACT
CATAGCCATCCTTGGTTCGGCTGGAATGACTTGCATCACCGTGTGTGGCGTTTCTGATTATGTTGCAA
CTGAAGAGAGCAAATGTCCAAGGAGAATGGCTCAGGCATTCCAGAACGTGAGAGAAGAACCAGCTGTGC
AGTTTAACTCAGGAACTCTGGCCCTAACAGGAAGGCCAAAAACAATCCGGATCCCACAATTTATCCTGT
GCTTGACTGGAATGACATCAAGTTTCAAGACGTGATCGGAGAGGGCAACTTTGGCCAGGTTCTGAAGGCA
CGCATCAAGAAGGATGGGTTACGGATGGATGCCGCCATCAAGAGGATGAAAGAGTATGCCTCCAAGATG
ATCACAGGGACTTCGCAGGAGAACTGGAGTTCTTTGTAAGTGGACACCATCCAACATCATTAACTCT
CTTGGGAGCATGTGAACCCGAGGCTATTTGTACCTAGCTATTGAGTATGCCCGCATGGAACCTCCTG
GACTTCTGCGTAAGAGCAGAGTGTAGAGACAGACCCTGCTTTTGCATCGCCAACAGTACAGCTTCCA
CACTGTCTCCCAACAGCTTCTTCAATTTGCTGCAGATGTGGCCCGGGGATGGACTACTTGAGCCAGAA
ACAGTTTATCCACAGGGACCTGGCTGCCAGAAACATTTAGTTGGTGAAGTACATAGCCAAATAGCA
GATTTTGGATTGTACGAGGTCAAGAAGTGTATGTGAAAAGACAATGGGAAGGCTCCCAGTGCCTTGGAT
TGGCAATCGAATCACTGAACTATAGTGTCTATAACAACCAACAGTGTCTGGTCTATGGTGTATTGCT
CTGGGAGATTGTTAGCTTAGGAGCACCCCTACTGCGGCATGACGTGCGCGGAGCTCTATGAGAAGCTA
CCCCAGGGCTACAGGCTGGAGAAGCCCTGAACTGTGATGATGAGGTGTATGATCTAATGAGACAGTGT
GGAGGGAGAAGCCTTATGAGAGACCATCAATTTGCCAGATATTGGTGTCTTAAACAGGATGCTGGAAAG
ACGGAAGACATACGTGAACACCACACTGTATGAGAAGTTACCTATGCAGGAATTGACTGCTCTGCGGAA
GAAGCAGCCTAG
    
```

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja1670_g01.zip

Restriction Sites: SgfI-MluI

ACCN: NM_013690

Insert Size: 3372 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_013690.2</u> , <u>NP_038718.2</u>
RefSeq Size:	4654 bp
RefSeq ORF:	3372 bp
Locus ID:	21687
Cytogenetics:	4 43.34 cM

Gene Summary:

Tyrosine-protein kinase that acts as cell-surface receptor for ANGPT1, ANGPT2 and ANGPT4 and regulates angiogenesis, endothelial cell survival, proliferation, migration, adhesion and cell spreading, reorganization of the actin cytoskeleton, but also maintenance of vascular quiescence. Has anti-inflammatory effects by preventing the leakage of proinflammatory plasma proteins and leukocytes from blood vessels. Required for normal angiogenesis and heart development during embryogenesis. Required for post-natal hematopoiesis. After birth, activates or inhibits angiogenesis, depending on the context. Inhibits angiogenesis and promotes vascular stability in quiescent vessels, where endothelial cells have tight contacts. In quiescent vessels, ANGPT1 oligomers recruit TEK to cell-cell contacts, forming complexes with TEK molecules from adjoining cells, and this leads to preferential activation of phosphatidylinositol 3-kinase and the AKT1 signaling cascades. In migrating endothelial cells that lack cell-cell adhesions, ANGPT1 recruits TEK to contacts with the extracellular matrix, leading to the formation of focal adhesion complexes, activation of PTK2/FAK and of the downstream kinases MAPK1/ERK2 and MAPK3/ERK1, and ultimately to the stimulation of sprouting angiogenesis. ANGPT1 signaling triggers receptor dimerization and autophosphorylation at specific tyrosine residues that then serve as binding sites for scaffold proteins and effectors. Signaling is modulated by ANGPT2 that has lower affinity for TEK, can promote TEK autophosphorylation in the absence of ANGPT1, but inhibits ANGPT1-mediated signaling by competing for the same binding site. Signaling is also modulated by formation of heterodimers with TIE1, and by proteolytic processing that gives rise to a soluble TEK extracellular domain. The soluble extracellular domain modulates signaling by functioning as decoy receptor for angiopoietins. TEK phosphorylates DOK2, GRB7, GRB14, PIK3R1, SHC1 and TIE1.[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1). Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.