

Product datasheet for MC223704

Ddb1 (NM_015735) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ddb1 (NM_015735) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Ddb1
Synonyms:	127kDa; AA408517; p127-Ddb1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>MC223704 representing NM_015735 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTCGTACAACACTACGTCGTAACGGCGCAGAAGCCTACTGCAGTGAACGGCTGTGTGACCGGACACTTTA
CTTCAGCGGAAGATCTGAACCTGTTGATTGCTAAAAACACAAGATTAGAGATCTATGTGGTCACTGCGGA
GGGTCTTCGGCCAGTCAAAGAGGTGGGAATGTACGGGAAGATTGCAGTCATGGAGCTCTTCAGGCCCAAG
GGGAGAGCAAAGACCTGCTGTTTCATCTTGACTGCTAAATACAATGCCTGTATCTTGGAGTATAAGCAGA
GTGGCGAAAGCATTGACATCATTACAAGAGCCCATGGCAATGTCCAGGACCGCATTGGTCCGCCCTCGGA
GACTGGCATCATTGGCATCATTGACCCTGAGTGCCGTATGATTGGCCTTCGACTCTACGATGGCCTTTTC
AAGGTCATTCCACTAGACCGAGACAATAAAGAGCTCAAGGCCCTTAACATCCGCCTAGAGGAACTGCACG
TCATTGATGTCAAGTTCCTGTATGGTTGTCAAGCACCTACCATTGCTTTGTCTACCAGGATCCTCAGGG
GCGGCATGTGAAAACCTATGAGGTGTCTCTTCGAGAGAAGGAGTTCAACAAGGCCCTTGAAACAGGAG
AATGTAGAAGCCGAAGCTTCAATGGTGATAGCAGTCCCAGAGCCTTTTGGAGGTGCCATCATCATTGGAC
AGGAATCAATCACCTATCATAACGGTGATAAATACCTGGCAATTGCTCCTCCTATCATTAAAGCAAAGCA
TATTGTGTGCCATAACCGGGTGGATCCTAATGGCTCGCGGTACCTTCTGGGAGATATGGAGGACCGGCTT
TTCATGCTGCTCTTAGAGAAGGAGGAGCAGATGGATGGCACTGTCAACCCTCAAGGACCTTCGAGTGGAA
TCCTCGGAGAGACCTCCATTGCTGAGTGCCCTAATCTTGTGATAATGGTGTGATTTGTGCGGCTCTCG
CCTAGGTGACTCCAGCTCGTAAAGCTCAATGTCGACAGCAACGAACAAGGCTCCTATGTAGTGGCCATG
GAAACCTTTACCAATTTAGGACCCATTGTGGATATGTGTGTGGTGGACCTGGAGAGGCAGGGACAGGGTC
AGCTGGTCACTTGTCTGGAGCTTCAAGGAAGGTTCTTGGGATCATCCGGAACGGCATTGGAATCCA
TGAGCATGCCAGCATTGACTTACCAGGCATCAAAGGTTTATGGCCTCTGCGCTCTGACCCGGTCTGGGAG
ACGGATGACACTTGGTGTCTCTTTTGTGGGCCAGACAAGGTTCTCATGCTAAACGGAGAGGAAGTGG
AAGAGACAGAAGCTGATGGCTTTGTGGACGATCAGCAGACTTTCTTTTGTGGCAAGCTGGCTCATCAGCA
ACTTATCCAGATCACCTCAGCATCTGTGAGGTTGGTGTCTCAAGAGCCCAAGCTTTGGTCAAGGAGTGG
AAAGAGCCTCAGGGCAAGAACATCAGTGTGGCCTCCTGTAAACAGCAGCCAGGTGGTTGTTGCCGTGGAA



[View online »](#)

```
GGGCACTATATTACCTTCAGATCCACCCTCAGGAGCTCCGGCAAATCAGCCACACAGAGATGGAACATGA
AGTGGCATGCTTGGATATCACCCATTAGGGGACAGCAATGGCCTCTCCCACTTTGTGCCATTGGACTT
TGGACAGACATCTCAGCTCGATCCTGAAGCTCCCATCTTTGAACTACTACACAAGGAGATGCTGGGTG
GAGAGATCATTCTCGATCCATCCTGATGACCCTTTGAGAGTAGCCACTACCTCTTTGTGCCTTGGG
AGATGGGGCTCTTTCTACTTTGGGCTCAACATCGAGACAGGTTACTGAGTGACCGTAAAAAGGTGACT
TTGGGCACCCAGCCACGGTGTGAGGACATCCGTTCTCTTTCTACCACCAACGTCCTTTGCTTGTCTG
ACCGGCCACTGTCATCTATAGCAGCAATCACAAGTTGGTCTTCTCCAATGTCAACCTCAAGGAGTGAA
CTATATGTGTCTCTCAACTCAGATGGCTATCCTGACAGTCTGGCATTGGCCAATAACAGCACTCTCACC
ATTGGCACCATCGATGAGATCCAGAAGCTCCATATTTCGCACAGTTCCCTCTATGAGTCTCCAGGAAGA
TCTGCTATCAGGAAGTGTCTCAGTGCTTTGGGGTCTTTCCAGCCGATTGAAGTCCAAGATAGCAGTGG
AGGCACTACTGCTTTGAGGCCAGTGCCAGCACCCAGGCCCTGTCCAGCAGCGTCAGCTCCAGAAACTT
TTCTCCAGCAGTACGGCTCCTCACGAGACTCCTTTGGAGAAGAAGTGGAGTGCACAACCTCCTCATCA
TTGACCAGCACACCTTTGAAGTGTTCACGCCATCAGTTTCTACAAAACGAGTATGCCCTCAGCCTGGT
CTCCTGCAAGTTGGGTAAGACCCCAACTTACTTCATTGTGGGCACAGCCATGGTATATCCTGAAGAG
GCAGAGCCCAAGCAGGTCGATTGTGGTCTTCCAGTATTCAGATGGAAAACACAGACTGTGGCTGAAA
AGGAAGTAAAAGGAGCTGTGACTCTATGGTGGAGTTAATGGGAAGCTGCTTGCCAGCATCAATAGCAC
GGTGCGGCTCTATGAGTGGACCACAGAAAAGGAGCTTCGCACTGAGTGTAACTACAACAACATCATG
GCCCTGTACCTGAAGACCAAGGGCGACTTCATCCTGGTGGGCGACCTCATGCGCTCAGTGTCTGCTTG
CCTACAAGCCTATGGAAGGAACTTTGAAGAGATTGCCGAGACTTTAATCCCAACTGGATGAGTGCAGT
GGAGATCTTGGATGATGACAATTTCTGGGAGCTGAAAATGCCTTTAACTTGTGTGTGCAGAAGGAC
AGTGCTGCCACCACTGATGAGGAGCGGCAGCACCTACAGGAGTTGGTCTCTCCACCTGGGCGAGTTTG
TCAATGTGTTCTGCCATGGCTCCCTAGTTATGCAGAATCTGGGCGAGGCTTCTACCCCTACACAGGGCTC
GGTGCTCTTTGGCACAGTCAACGGCATGATAGGGCTGGTGACTCGCTGTCAGAGAGCTGGTATAACCTC
CTGTTGGACATGCAGAACCGACTCAATAAGGTCAATAAGGTCATCAAAAAGCGTGGGCAAGATTGAGCACTCCTTGGGA
GATCCTTTACACTGAGCGAAAAGACAGAACCCAGCCACAGGCTTCATCGATGGTGTCTGATTGAAAGTTT
CCTAGATATCAGCCGCCCTAAGATGCAGGAGGTTGTGGCAAACCTTGCAGTATGATGATGGCAGTGGTATG
AAGCGGGAGGCAACTGCAGATGACCTCATCAAAGTCGTGGAGGAGCTAACTCGGATCCAT
```

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

- Restriction Sites:** Sgfl-MluI
- ACCN:** NM_015735
- Insert Size:** 3423 bp
- OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_015735.1](#), [NP_056550.1](#)

RefSeq Size: 4171 bp

RefSeq ORF: 3423 bp

Locus ID: 13194

UniProt ID: [Q3U1J4](#)

Cytogenetics: 19 6.66 cM

Gene Summary: Required for DNA repair. Binds to DDB2 to form the UV-damaged DNA-binding protein complex (the UV-DDB complex). The UV-DDB complex may recognize UV-induced DNA damage and recruit proteins of the nucleotide excision repair pathway (the NER pathway) to initiate DNA repair. The UV-DDB complex preferentially binds to cyclobutane pyrimidine dimers (CPD), 6-4 photoproducts (6-4 PP), apurinic sites and short mismatches. Also appears to function as a component of numerous distinct DCX (DDB1-CUL4-X-box) E3 ubiquitin-protein ligase complexes which mediate the ubiquitination and subsequent proteasomal degradation of target proteins. The functional specificity of the DCX E3 ubiquitin-protein ligase complex is determined by the variable substrate recognition component recruited by DDB1. DCX(DDB2) (also known as DDB1-CUL4-ROC1, CUL4-DDB-ROC1 and CUL4-DDB-RBX1) may ubiquitinate histone H2A, histone H3 and histone H4 at sites of UV-induced DNA damage. The ubiquitination of histones may facilitate their removal from the nucleosome and promote subsequent DNA repair. DCX(DDB2) also ubiquitinates XPC, which may enhance DNA-binding by XPC and promote NER. DCX(DTL) plays a role in PCNA-dependent polyubiquitination of CDT1 and MDM2-dependent ubiquitination of TP53 in response to radiation-induced DNA damage and during DNA replication. DCX(ERCC8) (the CSA complex) plays a role in transcription-coupled repair (TCR). May also play a role in ubiquitination of CDKN1B/p27kip when associated with CUL4 and SKP2 (By similarity). The DDB1-CUL4A-DTL E3 ligase complex regulates the circadian clock function by mediating the ubiquitination and degradation of CRY1 (PubMed:26431207). DDB1-mediated CRY1 degradation promotes FOXO1 protein stability and FOXO1-mediated gluconeogenesis in the liver (PubMed:28790135).[UniProtKB/Swiss-Prot Function]