

Product datasheet for MC225585

Smim22 (NM_001253802) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Smim22 (NM_001253802) Mouse Untagged Clone
Tag: Tag Free
Symbol: Smim22
Synonyms: EG432995; Gm5480
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC225585 representing NM_001253802
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGTAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGC

ATGGGGCTTTCTGATGAGGAGCTGAGGCAGGAGCTGCAGGCCACTGCAGAGGAAGTACTGGGGAGGCTGA
 GGAGCCGGCAGCTCTCCAGTCGGAATGGGACGTCGCTGCCTTTGTGGTCTTTCTCACATTTGTGGGCAC
 TGTGCTGTTGCTGTTGCTGGTCTTTGTCCATTGCTGCTGTTGTTGCTGCTGCAACACCTCCCCCAGG
 CCCCAGGAAGGTGAGCCCCGGAAGGAAAGCGCAATGGTGTGGATAACTGGCCCTGGAACCTTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001253802

Insert Size: 276 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001253802.1, NP_001240731.1</u>
RefSeq Size:	385 bp
RefSeq ORF:	276 bp
Locus ID:	432995
Cytogenetics:	16 A1