

## Product datasheet for MC225611

### Rpl37rt (NM\_001271591) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Rpl37rt (NM_001271591) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Rpl37rt
Synonyms:	EyeLinc4; Gm13826; OTTMUSG00000014624
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC225611 representing NM_001271591 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGCATCGCC

ATGACGAAGGGAACGTCATCCTTTGGTAAGCGTCGCAACAAGACGCACACGCTGTGCCGCCGCTGTGGCT  
 CCAAGGCTACCACTTCAGAAGTCGACTTGTGGCAAGTGTGGCTACCTGCCAAGCGCAAGAGGAAGTA  
 TAACTGGAGTGCCAAGGCTAAGAGACGAAACACTACCGGGACTGGTCGGATGAGGCACCTAAAGATTGTC  
 TACCGCAGATTGACATGGATTCCGTGAGGGAACAACACCGAAACCAAGAGGGCAGCTGTCGAGCAT  
 CCAGTTCTTCTTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001271591
Insert Size:	294 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.


[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>Note:</b>	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
<b>RefSeq:</b>	<u><a href="#">NM_001271591.1</a></u> , <u><a href="#">NP_001258520.1</a></u>
<b>RefSeq Size:</b>	677 bp
<b>RefSeq ORF:</b>	294 bp
<b>Locus ID:</b>	100502825
<b>UniProt ID:</b>	<u><a href="#">Q9D823</a></u>
<b>Cytogenetics:</b>	5 55.99 cM