

## Product datasheet for **MC225615**

### Hikeshi (NM\_001291289) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Hikeshi (NM\_001291289) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Hikeshi  
**Synonyms:** 0610007P06Rik; 1110002N09Rik; l(7)6Rn; l7Rn6  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC225615 representing NM\_001291289  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGAATATTGTGCGAACCCCATCTGTTGCCAGATTGGCATCTCGGTGGAATTGTTGGACAGTCTGGCTC  
AGCAGACTCCTGTAGGCAGTGTGCCGTGCTCCTCCGTTGACTCCTTTACGCAGTTCACACAGAAGATGTT  
GGACAACCTCTACAATTTTGCTTCATCATTGCTCTCTCTCAGGCCAGATGACACCAAATCCATCAGAA  
ATGTTTCATCCCAGCAAATGTGGTTCTGAAATGGTATGAAAACCTTCAAAGACGACTAGCACAGAACCCCC  
TCTTTTGAAAACATAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** Sgfl-MluI  
**ACCN:** NM\_001291289  
**Insert Size:** 297 bp  
**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).  
**OTI Annotation:** Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.



<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001291289.1, NP_001278218.1</u>
<b>RefSeq Size:</b>	847 bp
<b>RefSeq ORF:</b>	297 bp
<b>Locus ID:</b>	67669
<b>Cytogenetics:</b>	7 50.4 cM
<b>Gene Summary:</b>	<p>Acts as a specific nuclear import carrier for HSP70 proteins following heat-shock stress: acts by mediating the nucleoporin-dependent translocation of ATP-bound HSP70 proteins into the nucleus. HSP70 proteins import is required to protect cells from heat shock damages. Does not translocate ADP-bound HSP70 proteins into the nucleus (By similarity). May also be indirectly required for organization and/or function of the secretory apparatus in Clara cells in lung.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (5) lacks an alternate exon and initiates translation at an alternate start codon compared to variant 1. The resulting isoform (d) is shorter at the N-terminus compared to isoform a. Variants 4 and 5 both encode the same isoform (d).</p>