

Product datasheet for **MC225637**

3830406C13Rik (NM_001284386) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	3830406C13Rik (NM_001284386) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	3830406C13Rik
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>MC225637 representing NM_001284386 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGTTACTTCAGAAAATGAAAAACAAATTTGGAGATGAAAACACAGAAAGGGCATCTCTCCTCCAAGCCA
CTGAGACTGCTTCTAGAAGAACTCCGACTCTTGAAGGACATAGATGCTGCAGAGGAGGCTTTTCAGAC
CAAGCTTATCCCACACCCACAGCCTTCGATGCTTCTCTTGAGACTCGTTACTGGGCATCAGTAGAAGAA
CATATCCCAAGTGGAAGTCTTTCTTTAGGAAGAGCACCATATCCAATTGGTGCTGAAAATCAAATG
AGGTACCATTGTTCAAACCGAGGCAAAGTTG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001284386
Insert Size:	315 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001284386.1, NP_001271315.1

RefSeq Size: 2338 bp

RefSeq ORF: 315 bp

Locus ID: 218734

UniProt ID: Q8BGD0

Cytogenetics: 14 A1