

Product datasheet for **MC225667**

Plac9b (NM_001270503) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Plac9b (NM_001270503) Mouse Untagged Clone
Tag: Tag Free
Symbol: Plac9b
Synonyms: ENSMUSG00000072674; Gm10393
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC225667 representing NM_001270503
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCAGGCGCTGCTCTGCGCGCTAGCCGGGCTCGCCCTGCTCCGTGCCGGGACGGGCGAATGGGGCCAAG
GTCCCAGGGACACCCCGGACGGCGAGCTGCCGAGTCCCCAGTCCCCCTGGAGATCTAGCTGGAGGCC
AGGCTGTGACAGACACGCGGCTGTGCAAAGCGGTTAGACATTATGGAGGAGACGGTGGAGAAGACAGTG
GAGCACCTGGAGGCGGAAGTGACAGGTCTGCTGGGCTGCTGGAGGAAGTGGCTCAAACCTCCCACAG
GGCCCTTCAGCCCAAACCTGACTTGCTTGAGATGATGGTTTCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001270503

Insert Size: 327 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.



[View online »](#)

| | |
|-------------------------------|---|
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001270503.1</u> , <u>NP_001257432.1</u> |
| RefSeq Size: | 1020 bp |
| RefSeq ORF: | 327 bp |
| Locus ID: | 100039246 |
| UniProt ID: | <u>Q8K262</u> |
| Cytogenetics: | 14 |