

Product datasheet for **MC225829**

Ly6a (NM_001271417) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ly6a (NM_001271417) Mouse Untagged Clone
Tag: Tag Free
Symbol: Ly6a
Synonyms: Ly-6A.2; Ly-6A/E; Ly-6E.1; Sca-1; Sca1; TAP
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC225829 representing NM_001271417
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGGACACTTCTCACACTACAAAGTCTGTTTGCTGATTCTTCTGTGGCCCTACTGTGTGCAGAAAGAG
CTCAGGGACTGGAGTGTTACCACTGCTATGGAGTCCCATTTGAGACTTCTTGCCCATCAATTACCTGCC
CTACCTGATGGAGTCTGTGTTACTCAGGAGGCAGCAGTTATTGTGGATTCTCAAACAAGGAAAGTAAAG
AACAATCTTTGCTTACCCATCTGCCCTCCTAATATTGAAAGTATGGAGATCCTGGGTACTAAGGTCAACG
TGAAGACTTCCTGTTGCCAGGAAGACCTCTGCAATGTAGCAGTCCCAATGGAGGCAGCACCTGGACCAT
GGCAGGGGTGCTTCTGTTTCAGCCTGAGCTCAGTCTCTGCAGACCTTGCT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001271417

Insert Size: 405 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001271417.1, NP_001258346.1</u>
RefSeq Size:	987 bp
RefSeq ORF:	405 bp
Locus ID:	110454
UniProt ID:	<u>P05533</u>
Cytogenetics:	15 34.29 cM
Gene Summary:	<p>T-cell activation.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) uses an alternate splice site in the 5' UTR, compared to variant 1. Variants 1, 2, 3, 4, 5, and 6 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>