

Product datasheet for **MC225838**

Plekhb1 (NM_001291336) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Plekhb1 (NM_001291336) Mouse Untagged Clone
Tag: Tag Free
Symbol: Plekhb1
Synonyms: AI315628; evt-1; PHR1; Phret1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC225838 representing NM_001291336
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGC**C

ATGGCCCTGGTGAGGGGGCGGCTGGCTGTGGAGACAGAGCTCCATCCTCCGCCGCTGGAAGCGGAATTGGT
 TCGCTGTGGTTGGATGGCACACTGGGTTACTACCATGACGAGACGGCACAGGACGAGGAGACCGAGT
 AGTCATCCACTTCAATGTCCGAGACATAAAGGTCGGCCAAGAGTGTGAGGATGTGCAGCCCCAGAGGGC
 AGGAGCCGAGATGGCCTGCTGACAGTGAACCTGCGGGAGGGTTCGCCCTGCACCTGTGCGCTGAGACCC
 GGGATGATGCCATGTACGCGTCTACAGCCCGTATCAGGACTACTATGAGGTGGTCCACCCAACGCACAT
 GAAGCCACATATGTCCGAGCTACTATGGGCCACCTTATGCAGGTCTGGTGTGA

AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-RsrII

ACCN: NM_001291336

Insert Size: 405 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001291336.1, NP_001278265.1</u>
RefSeq Size:	1784 bp
RefSeq ORF:	405 bp
Locus ID:	27276
UniProt ID:	<u>Q9QYE9</u>
Cytogenetics:	7 E2