

Product datasheet for MC225866

9330182014Rik (NM_001256056) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	9330182014Rik (NM_001256056) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	9330182014Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC225866 representing NM_001256056

Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGTATCAGAGAAATGGGGAGTTACTTCCGGTTGTTTGGACCTCCCAAATTATACTCTCCAGAATTTA
 CAGAGAAAGCCAGACAGTGTCTTCACTGCCTCCTGACTGCTAGAATGTCTACCCTGTGTGCTTACCATGA
 AGTAAAGAAACCCCTCTCCATTATACATTATTCCTGTCAACATGATGTTCTGCCAAGGATATGAGACCA
 AGCAACACAGACAGGATCCTCTCATATGAACTGTTAAGACGATGGAGACTCCTGGAGATGGACTCAATGC
 ATTTGCTTCACGAGGAGTCTTTGGGATGTTATGGGATATTATGGTTGACCAGGGCATGTTCCACAGACT
 CAGGTATTGGTGCTTTAGCTGGCAGCTTCAGGGGTTTTTCAGAAGTGATTGGATCCTAAGAAGCTGATG
 TAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001256056
Insert Size:	423 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001256056.1</u> , <u>NP_001242985.1</u>
RefSeq Size:	2053 bp
RefSeq ORF:	423 bp
Locus ID:	328531
Cytogenetics:	15 B3.1